

Lecture 4

Registry of HLA epitopes and
HLAMatchmaker analysis for
antibody-verified epitopes

16th Workshop Project: International HLA Epitope Registry

GOALS:

- Develop a notation system for epitopes encoded by HLA-ABC, DRB1/3/4/5, DQAB, DPAB and MICA
- Establish a database of antibody-verified HLA epitopes
- Offer clinically relevant search functions

The registry does not address so-called cellular HLA epitopes defined by alloreactive T-cells

The Website-Based International Registry of Antibody-Verified HLA Epitopes

Int J. Immunogenetics, 40: 54-59

<http://www.epregistry.com>

Developed and managed by

Marilyn Marrari and Rene Duquesnoy

University of Pittsburgh School of Medicine, Pittsburgh PA

and

Luiz Claudio Sousa, Adalberto da Silva, Keylla Maria
Aita, Jose Renato Barroso, Mário Sérgio Marroquim
and Semiramis do Monte

Federal University of Piaui, Teresina, PI, Brazil

HLA Epitope Registry: www.epregistry.ufpi.br

(can be viewed in Microsoft Internet Explorer,
Google Chrome or Mozilla Firefox browsers)

- Website information is for research use only. Users must agree to terms of use before using website:



WEBSITE TERMS OF USE

WEBSITE DISCLAIMER

By using the site, you agree to the following terms of use, which constitute an agreement between the website and you.

These terms may be changed from time to time. You may view the most current version through the link that appears on the bottom of all pages on the website.

1. This website is designed for research purposes only. The contents are not intended for making clinical decisions regarding donor selection or patient care.
2. The website disclaims any liability whatsoever for any documentation, information, programs, software or other material that is or may become part of the website. We do not warrant or guarantee that the information will not be offensive to any user and/or may not meet your needs and requirements.
3. The website may include links to third party Web service providers. These links are provided solely as a convenience to you, and not as an endorsement of the content on such sites. We are not responsible for the content of linked third-party sites.

I have read the Terms of Use shown above and agree with them.

I agree

I Disagree

Registry Home Page

The screenshot shows the HLA Epitope Registry website. The browser address bar displays 'epregistry.ufpi.br'. The website has a blue and green abstract banner at the top. Below the banner is a navigation bar with links: Home (active), Databases, Resources, and Contact Us. On the left side, there is a sidebar with a 'Description' link and a list of database links: Registration, New Epitopes, ABC Database, DRB Database, DQ Database, DP Database, and MICA Database. The main content area features a 'Description' section with a heading, a paragraph about the website's research purpose, a list of five databases, and a numbered list of seven layout features. At the bottom, a paragraph explains the data organization and provides a link for more information.

epregistry.ufpi.br

HLA Epitope Registry

Home Databases Resources Contact Us

Description

This website is designed for research purposes only. The contents are not intended for making clinical decisions regarding donor selection or patient care.

There are five separate databases: ABC, DRB1/3/4/5, DQB + DQA, DPB + DPA and MICA.

Their layouts display:

- 1 Epitope names
- 2 Polymorphic residue descriptions
- 3 Epitope frequencies
- 4 Antibody reactivity descriptions of "confirmed" or "provisional" antibody-defined epitopes
- 5 Information about corresponding "structural" epitopes
- 6 Epitope-carrying alleles in Luminex panels
- 7 Listings of all alleles with antibody-defined epitopes.

Each epitope has its own row and epitopes are sorted according their sequence positions. Under each epitope there are rows for possible variants with distinct molecular configurations. Each epitope database has search options to identify repertoires of antibody-defined epitopes on selected alleles and epitopes that are mismatched for a given HLA type. For more detailed information about the HLA Epitope Registry, [click here](#).

Class I Eplet Examples

| Eplet | Description | Eplet-Carrying Alleles in Luminex Panel |
|-------------------|--|--|
| 62GE | 62G,63E,65R | A*02:01/02/03/05/06,B*57:01/03,*58:01 |
| 166DG | 166D,167G | A*01:01,*23:01,*24:02,*80:01,B*15:12 |
| 163EW | 162G163E166E167W | A*66:02,B*07:02*13:01/02,*27:05/08,*40:01/02/06,*47:01,*48:01,*73:01,*81:01,C*02:02,*17:01 |
| 44KM ₃ | 43Q44K45M46E / 149A150V151H152A / 158V | A*01:01, A*36:01 |

| Locus | Eplet | Description | Eplet-Carrying Alleles in Luminex Panel |
|-------|-------|--------------|---|
| MICA | 24T | 24T, 26V,36C | A*001,*012,*018 |
| MICA | 175S | 173E,175S | A*004,*006,*009,*016,*019,*033,*036 |

Class I ABC Database

Databases

ABC

DRB

DQB + DQA



DPB + DPA

MICA


► [Filters](#)

| Epitope | Polymorphic Residues | Frequency | Antibody Reactivity | Struct Epitope | Luminex Alleles | All Alleles |
|---------|----------------------|-----------|---------------------|----------------------|--|--|
| 163LE | 163L166E | | | | B*15:01, B*15:02, B*15:03, B*15:10, B*15... More | A*26:04, B*07:20, B*07:24, B*07:60, B*08... More |
| 163R | 163R | | Provisional | View | A*01:01, A*11:01, A*11:02, A*25:01, A*26... More | A*01:01, A*01:02, A*01:03, A*01:06, A*01... More |
| 163RD | 162G163R166D167G | | Confirmed | View | A*01:01 More | A*01:01, A*01:02, A*01:03, A*01:06, A*01... More |
| 163RE | 162G163R166E167W | | Confirmed | View | A*11:01, A*11:02, A*25:01, A*26:01, A*43... More | A*02:101, A*02:135, A*02:154, A*02:38, A... |

Antibody Reactivity Report for 163RE

| Epitope: | 163RE (confirmed) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|-----------------------------------|---|--------|--------|--------|--------|------------------|----|----|----|------------------|---|---|---|----------------|-------|------|--|----------------|-------|------|--|----------------|-------|------|--|-----------------------------------|-------|------|--|----------------|------|------|--|----------------|-------|-------|--|--------------|---------|-------|--|------------------------------|----------|-------|--|---------|-----------------|-----|-------|---|---|-------|---|---|-------|---|---|-------|---|---|-------|---|---|-------|---|---|-------|---|---|-------|---|---|-------|---|---|-------|---|---|
| Luminex single allele reactivity: | <p>Human monoclonal antibody VP5G3 (IgM) from Arend Mulder (Leiden): Sensitization due to pregnancy. Antibody producer types as HLA-A24, A32; B18, B60; Cw10, Cw7. Immunizing haplotype: HLA-A26; B55 (Bw6); Cw1 Data supplied by Marilyn Marrari (Pittsburgh):</p> <table><thead><tr><th></th><th>PANEL</th><th>OL MFI</th><th>GP MFI</th></tr></thead><tbody><tr><td>Positive control</td><td>nt</td><td>nt</td><td>nt</td></tr><tr><td>Negative control</td><td>9</td><td>5</td><td>5</td></tr><tr><td>163RE+ A*11:01</td><td>11150</td><td>8163</td><td></td></tr><tr><td>163RE+ A*11:02</td><td>12567</td><td>4246</td><td></td></tr><tr><td>163RE+ A*25:01</td><td>13029</td><td>7194</td><td></td></tr><tr><td>Immunizing allele: 163RE+ A*26:01</td><td>12521</td><td>9514</td><td></td></tr><tr><td>163RE+ A*43:01</td><td>7510</td><td>1165</td><td></td></tr><tr><td>163RE+ A*66:01</td><td>12498</td><td>11044</td><td></td></tr><tr><td>Self Alleles</td><td>21 ± 11</td><td>9 ± 4</td><td></td></tr><tr><td>Other 163RE-negative alleles</td><td>84 ± 278</td><td>9 ± 3</td><td></td></tr></tbody></table> <p>Data from El-Awar et al, <i>Human Immunology</i> 68: 170-180, 2007: Serum eluate Z2076.@0 (absorbed with A*25:01 cells), no information on HLA types of antibody producer/immunizer</p> <p>Epitope #209 Z2076.@0: Eluted from A2501 rHLA cell line</p>  <table><thead><tr><th>Antigen</th><th>aa position 163</th><th>166</th></tr></thead><tbody><tr><td>A0101</td><td>R</td><td>D</td></tr><tr><td>A1101</td><td>R</td><td>E</td></tr><tr><td>A1102</td><td>R</td><td>E</td></tr><tr><td>A2501</td><td>R</td><td>E</td></tr><tr><td>A2601</td><td>R</td><td>E</td></tr><tr><td>A4301</td><td>R</td><td>E</td></tr><tr><td>A6601</td><td>R</td><td>E</td></tr><tr><td>B4201</td><td>T</td><td>E</td></tr><tr><td>A3401</td><td>T</td><td>E</td></tr><tr><td>A6602</td><td>E</td><td>E</td></tr></tbody></table>  | | PANEL | OL MFI | GP MFI | Positive control | nt | nt | nt | Negative control | 9 | 5 | 5 | 163RE+ A*11:01 | 11150 | 8163 | | 163RE+ A*11:02 | 12567 | 4246 | | 163RE+ A*25:01 | 13029 | 7194 | | Immunizing allele: 163RE+ A*26:01 | 12521 | 9514 | | 163RE+ A*43:01 | 7510 | 1165 | | 163RE+ A*66:01 | 12498 | 11044 | | Self Alleles | 21 ± 11 | 9 ± 4 | | Other 163RE-negative alleles | 84 ± 278 | 9 ± 3 | | Antigen | aa position 163 | 166 | A0101 | R | D | A1101 | R | E | A1102 | R | E | A2501 | R | E | A2601 | R | E | A4301 | R | E | A6601 | R | E | B4201 | T | E | A3401 | T | E | A6602 | E | E |
| | PANEL | OL MFI | GP MFI | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Positive control | nt | nt | nt | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Negative control | 9 | 5 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 163RE+ A*11:01 | 11150 | 8163 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 163RE+ A*11:02 | 12567 | 4246 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 163RE+ A*25:01 | 13029 | 7194 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Immunizing allele: 163RE+ A*26:01 | 12521 | 9514 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 163RE+ A*43:01 | 7510 | 1165 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 163RE+ A*66:01 | 12498 | 11044 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Self Alleles | 21 ± 11 | 9 ± 4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Other 163RE-negative alleles | 84 ± 278 | 9 ± 3 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Antigen | aa position 163 | 166 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A0101 | R | D | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A1101 | R | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A1102 | R | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A2501 | R | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A2601 | R | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A4301 | R | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A6601 | R | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| B4201 | T | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A3401 | T | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A6602 | E | E | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Other antibody assays: | <p>Data from Mulder et al, <i>Tissue Antigens</i> 52: 393-396, 1998: VP5G3 reacts with HLA-A11, A25, A26, A66 in CDC assays; A43 not tested.</p> <p>Data from Fernandez-Vina et al, in <i>Immunobiology of the Human MHC vol. I</i>, John A. Hansen (Ed), pp 890-931, 2006: VP5G3 (13-WS 0061) positive with HLA-A*11:01, *11:02, *11:03, *24:10, *25:01, *25:02, *26:01, *26:02, *26:03, *26:08, *26:15, *66:01 in CDC assays in 13th IHWS. No HLA-A43 cells were tested.</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Comments: | Equivalent to TerEp #209. Epitope is not shared by HLA-A*66:02, but is shared by *24:10, which was positive in CDC assays. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Epitope-sharing Luminex alleles


 Epitope Registry

[Home](#) [Database](#)

| | | | | | |
|--------|------------------|--|-----------|--|--|
| 163DLS | 162D163L166 | | | | |
| 163E | 163E | | | | 0:0 7:0 |
| 163EW | 162G163E166E167W | | Confirmed | | B*13... More A*66:02, B*07:0 B*07:03, B*13:0 B*13... More |

Luminex Alleles of Epitope: 163RE Confirmed ×
A*11:01, A*11:02, A*25:01, A*26:01, A*43:01, A*66:01
[Close](#)

All Alleles that Share Epitope

 Epitope Registry

[Home](#) [Database](#)

| | |
|--------|-------------|
| 163EW | 162G163E166 |
| 163GLS | 162G163L166 |
| 163LW | 162G163L166 |

Alleles of Epitope: 163RE Confirmed ×

A*02:101, A*02:135, A*02:154, A*02:38, A*11:01, A*11:02, A*11:03, A*11:05, A*11:06, A*11:07, A*11:08, A*11:09, A*11:10, A*11:11, A*11:12, A*11:13, A*11:14, A*11:15, A*11:16, A*11:17, A*11:18, A*11:19, A*11:20, A*11:22, A*11:23, A*11:24, A*11:25, A*11:26, A*11:29, A*11:30, A*11:31, A*11:32, A*11:33, A*11:34, A*11:36, A*11:37, A*11:40, A*11:41, A*11:42, A*24:10, A*25:01, A*25:02, A*25:03, A*25:04, A*25:05, A*25:07, A*25:08, A*26:01, A*26:02, A*26:03, A*26:05, A*26:06, A*26:07, A*26:08, A*26:10, A*26:12, A*26:13, A*26:14, A*26:15, A*26:16, A*26:17, A*26:18, A*26:19, A*26:20, A*26:21, A*26:22, A*26:23, A*26:24, A*26:26, A*26:27, A*26:28, A*26:30, A*26:31, A*26:32, A*26:33, A*26:35, A*26:36, A*26:37, A*26:38, A*31:24, A*36:04, A*43:01, A*66:01, A*66:04, A*66:05, A*66:06, A*66:07, A*66:08, A*66:09

Close

Select a Database to View - DRB

Databases

ABC

DRB

DQB + DQA

DPB + DPA

MICA

› Filters

| Epitope | Polymorphic Residues | Frequency | Antibody Reactivity | Struct Epitope | Luminex Alleles | All Alleles |
|---------|----------------------|-----------|---------------------|----------------|--|--|
| 4Q | 4Q | | Provisional | | DRB1*07:01, DRB1*09:01, DRB4*01:01, D... More | DRB1*07:01, DRB1*07:09, DRB1*09:01, D... More |
| 4R | 4R | | | | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More |
| 6C2 | 6C(157I) | | | | DRB5*02:02 More | DRB1*13:112, DRB1*13:128, DRB5*02:02,... More |
| 6R2 | 6R(157T) | | | | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More |
| 12TSE | 12T13S14E | | Confirmed | | DRB1*03:01, DRB1*03:02, DRB1*03:03, D... More | DRB1*03:01, DRB1*03:02, DRB1*03:03, D... More |
| 13FE | 13F14E | | Provisional | | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More |
| 13FEL | 13F14E26L | | Provisional | | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More | DRB1*01:01, DRB1*01:02, DRB1*01:03, D... More |

Epitope Database Search Options

1. Identify the epitope repertoire of any selected allele (Luminex or non-Luminex)

Choose antibody-verified or all epitopes

2. Identify all mismatched epitopes for a given HLA phenotype

Choose antibody-verified or all epitopes

These options can be combined to identify a donor allele's repertoire of epitopes that are mismatched for a recipient's HLA type

Identifying epitopes on non-Luminex alleles mismatched for a patient

Enter the non-Luminex allele of interest (e.g., donor) in the “Search by allele” field and the patient HLA type in the “HLA type of recipient” field; results show one epitope on non-Luminex A*03:02 that is not shared by the common A*03:01 found in the Luminex panels.

Filters

Filters parameters: Antibody Reactivity is All and HLA Typing is **A*03:02** and HLA Type of recipient is A*01:01, A*02:03, B*07:02, B*27:05, C*02:02, C*07:02

Back

non-Luminex allele on potential donor

epitope shared by non-Luminex A*03:02 but not common A*03:01

| Epitope | Polymorphic Residues | Frequency | Antibody Reactivity | Struct Epitope | Luminex Alleles | All Alleles |
|-------------|----------------------|-----------|---------------------|----------------|--|--|
| 66NAQ | 66N69A70Q | | | | A*03:01, A*11:01, A*11:02, A*29:01, A... More | A*01:13, A*01:17, A*02:103, A*02:56... More |
| 66NV | 66N67V69A | | | | A*03:01, A*11:01, A*11:02, A*25:01, A... More | A*01:07, A*01:23, A*02:08, A*02:103... More |
| 69AQS | 69A70Q71S | | | | A*03:01, A*11:01, A*11:02, A*29:01, A... More | A*01:13, A*01:17, A*02:34, A*02:35, A... More |
| 70QS | 70Q71S | | | | A*03:01, A*11:01, A*11:02, A*29:01, A... More | A*01:13, A*01:17, A*02:103, A*02:34... More |
| 73TDVD | 73T74D76V77D | | | | A*03:01, A*11:01, A*11:02, A*30:01, A... More | A*01:13, A*01:28, A*02:35, A*02:48, A... More |
| 149AAH | 149A150A151H | | Confirmed | | A*02:01, A*02:02, A*02:05, A*02:06, A... More | A*01:12, A*01:19, A*01:21, A*01:26, A... More |
| 151HV | 150A151H152V | | | | A*02:01, A*02:02, A*02:05, A*02:06, A... More | A*01:12, A*01:19, A*01:21, A*02:01, A... More |
| 156Q | 156Q | | | | A*11:01, A*11:02, A*24:02, A*24:03, C... More | A*01:12, A*01:19, A*01:25, A*02:12, A... More |

If screening results show antibody reactivity towards epitope 156Q, A*03:02 is likely an unacceptable allele!

Combining Options

Use Filters 2 and 3 in combination to display epitopes on an allele that are mismatched for a given HLA type, e.g., to identify epitope mismatches for a patient's HLA phenotype that are found on a specific donor allele.

Optional filter to refine epitope search:

Search by allele:

A*11:01

Example: A*01:01

HLA Type of recipient: If an allele is unknown, leave the space blank.

| A | A | B | B | C | C |
|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|
| A*02:01 | A*03:01 | B*07:02 | B*35:01 | C*07:02 | |
| Example: A*01:01 | Example: A*29:02 | Example: B*08:01 | Example: B*40:01 | Example: C*03:04 | Example: C*07:01 |

Q Filter

Press 'Filter' button to see results

Enter single allele and HLA phenotype information

Results of combined options

Filters

Filters parameters: Antibody Reactivity is All and HLA Typing is A*11:01 and HLA Type of recipient is A*02:01, A*03:01, B*07:02, B*35:01, C*07:02

Back

| Epitope | Polymorphic Residues | Frequency | Antibody Reactivity | Struct Epitope | Luminex Alleles | All Alleles |
|---------|----------------------|-----------|---------------------|----------------------|---|---|
| 151AHA | 150A151H152A | | Confirmed | | A*11:01, A*11:02 More | A*01:26, A*02:27, A*11:01, A*11:02, A... More |
| 151HA | 151H152A | | | | A*01:01, A*11:01, A*11:02, A*36:01 More | A*01:01, A*01:02, A*01:03, A*01:06, A... More |
| 156Q | 156Q | | | | A*11:01, A*11:02, A*24:02, A*24:03, C... More | A*01:12, A*01:19, A*01:25, A*02:12, A... More |
| 163R | 163R | | Provisional | | A*01:01, A*11:01, A*11:02, A*25:01, A... More | A*01:01, A*01:02, A*01:03, A*01:06, A... More |
| 163RE | 162G163R166E167W | | Confirmed | View | A*11:01, A*11:02, A*25:01, A*26:01, A... More | A*02:101, A*02:135, A*02:154, A*02:38... More |

We found 5 epitopes in the database.

The HLA Epitope Registry is a Work-In-Progress

- As of July 2013, 85 HLA-ABC, 54 DRB1/3/4/5, 21 DQ (18 DQB1, 3 DQA1), 7 DPB1 and 24 MICA antibody-verified epitopes are listed
- Each database will be updated on a quarterly basis
- Currently working on a tool to allow users to print or create and download pdf files of results

New Antibody-Verified Epitopes

- An important aspect is the participation by HLA professionals who have identified antibodies specific for new and not-so-well described epitopes
- From the website one can download instructions for how to submit information about epitope-specific antibody reactivity

HLAMatchmaker Analysis of Sera from Sensitized Patients

Use of antibody-verified class I epitopes
recorded in the International HLA Epitope
Registry

Analysis of Epitope Specificity of HLA Antibody

HLA information about the sensitizing event

- High-resolution HLA typing of patient and the immunizing HLA antigen
- Provides information about patient's exposure to mismatched eplets
- Identification of immunizer eplets on alleles in the Luminex panel

Everybody Likes the Luminex
Assay with Single Alleles

but

It's Not a Perfect Test!

Potential Limitations of Single Allele Luminex Assays

- Some Luminex kits have certain alleles with significant MFI values unrelated to allosensitization events (technical problems or do they reflect “natural” antibodies?)

Unexplained high MFI values with 129 first pregnancy sera with no antibodies against paternal HLA class I alleles (Study with Stefan Schaub, Basel, Switzerland)

| | <u>Nr of cases</u> | <u>Mean MFI + SD (Range)</u> | | | |
|------------------|--------------------|------------------------------|--|--|--|
| Self alleles | 129 | 6 ± 7 (1-55) | | | |
| Paternal alleles | 299 | 9 ± 17 (1-108) | | | |

| <u>Luminex allele</u> | <u>Nr of cases</u> | <u>Adjusted MFI*</u> | <u>Luminex allele</u> | <u>Nr of cases</u> | <u>Adjusted MFI*</u> |
|-----------------------|--------------------|----------------------------|-----------------------|--------------------|----------------------|
| A*11:02 | 6 | 221,254,428,1033,1122,1493 | B*14:02 | 1 | 835 |
| B*15:12 | 5 | 201,209,314,391,709 | B*40:02 | 1 | 365 |
| A*29:01 | 4 | 237,251,519,3874 | B*42:01 | 1 | 333 |
| A*29:02 | 4 | 204,236,1864,3819 | B*15:16 | 1 | 317 |
| B*45:01 | 3 | 305,526,1718 | B*67:01 | 1 | 310 |
| A*68:02 | 3 | 304,889,1141 | B*13:02 | 1 | 300 |
| B*73:01 | 3 | 177,320,629 | B*44:03 | 1 | 290 |
| B*58:01 | 2 | 222,342 | B*54:01 | 1 | 258 |
| A*80:01 | 2 | 880,1726 | B*18:01 | 1 | 225 |
| B*57:01 | 2 | 203,229 | B*44:02 | 1 | 225 |
| B*82:01 | 2 | 272,578 | A*26:01 | 1 | 224 |
| C*17:01 | 1 | 2113 | C*07:02 | 1 | 211 |
| C*03:02 | 1 | 922 | A*25:01 | 1 | 203 |

* Adjusted MFI = observed MFI - (Mean MFI with self + 3 StDev) Only >200 values are shown

Antibody Verification of HLA Class I Epitopes

- Duquesnoy RJ, Marrari M, Mulder A, da M. Sousa LCD, Da Silva AS, do Monte SJH. 'First Report on the Antibody Verification of HLA-ABC Epitopes Recorded in the HLA Epitope Registry'. *Tissue Antigens* 83:391-400, 2014
 - 97 HLA-ABC antibody-verified epitopes have been recorded, 62 correspond to eplets and 35 are defined by eplets paired with other residue configurations

Serum Analysis with HLAMatchmaker

Three types of epitopes in the International Registry:

1. Antibody-verified epitopes defined by eplets*
2. Antibody-verified epitopes defined by eplet pairs*
3. Theoretical epitopes with distinct residue structures (might become antibody-verified?)

* Duquesnoy RJ, Marrari M, Mulder A, da M. Sousa LCD, Da Silva AS, do Monte SJH. 'First Report on the Antibody Verification of HLA-ABC Epitopes Recorded in the HLA Epitope Registry'. *Tissue Antigens* 83:391-400, 2014

97 HLA-ABC antibody-verified epitopes have been recorded, 62 correspond to eplets and 35 are defined by eplets paired with other residue configurations

HLAMatchmaker Analysis for Antibody-Verified Epitopes

- Class I SAB data on 28 first pregnancy sera provided by Stefan Schaub (Basel, Switzerland)
- High-resolution types of mother and paternal haplotype of child
- Analysis only done for sera whereby at least one immunizing allele had MFI>2500

| Case | HLA Type of Antibody Producer | | | MFI self | MFI Cut-off | MFI Pos Cont | Mismatched Immunizer Alleles | | | Number of Reactive Alleles | | | cPRA |
|------|-------------------------------|---------|---------|----------|-------------|--------------|------------------------------|--------------|--------------|----------------------------|-----------|------------|------|
| | | | | | | | MFI Values | | | HLA-A | HLA-B | HLA-C | |
| | | | | | | | | | | MFI Range | MFI Range | MFI Range | |
| #186 | A*11:01 | A*24:02 | B*07:02 | | | | A*03:02 | B*27:05 | C*01:02 | 1 | 7 | 2 | 32% |
| #186 | B*39:06 | C*07:02 | C*07:02 | 7±5 | 100 | 7002 | nd (01=286) | 8953 | 7245 | 286 | 104-8953 | 691-7245 | |
| #254 | A*01:01 | A*03:01 | B*40:02 | | | | A*31:01 | B*27:05 | C*01:02 | 0 | 2 | 8 | 56% |
| #254 | B*57:01 | C*02:02 | C*06:02 | 4±3 | 100 | 8855 | 2 | 11 | 3025 | | 2182-3061 | 210-3202 | |
| #247 | A*03:01 | A*03:01 | B*15:01 | | | | A*02:01 | Match | C*03:03 | 13 | 3 | 0 | 71% |
| #247 | B*35:01 | C*04:01 | C*07:02 | 15±29 | 200 | 9736 | 4347 | | 1 | 314-5951 | 724-1042 | | |
| #233 | A*02:01 | A*23:01 | B*15:18 | | | | Match | B*40:02 | C*02:02 | 2 | 31 | 1 | 71% |
| #233 | B*51:01 | C*07:04 | C*15:02 | 86±39 | 200 | 8285 | | 9647 | 430 | 7976 | 226-11608 | 430 | |
| #373 | A*02:01 | A*11:01 | B*35:03 | | | | A*29:01 | B*08:01 | C*07:02 | 25 | 32 | 3 | 73% |
| #373 | B*35:01 | C*04:01 | C*04:01 | 1±0 | 100 | 8695 | 8649 | 10629 | 5 | 559-10276 | 117-10629 | 112-144 | |
| #129 | A*01:01 | A*23:01 | B*15:17 | | | | A*30:02 | B*13:02 | C*06:02 | 28 | 27 | 7 | 75% |
| #129 | B*49:01 | C*07:01 | C*07:01 | 3±3 | 100 | 9719 | 10888 | 8982 | 5712 | 118-11108 | 103-8982 | 256-7653 | |
| #355 | A*32:01 | A*32:01 | B*15:01 | | | | A*11:01 | B*51:01 | C*15:02 | 1 | 13 | 9 | 76% |
| #355 | B*38:01 | C*03:03 | C*12:03 | 46±45 | 200 | 7861 | 47 | 8272 | 11977 | 583 | 243-8272 | 1075-11977 | |
| #354 | A*01:01 | A*02:01 | B*35:03 | | | | A*30:01 | B*51:01 | C*14:03 | 8 | 20 | 0 | 77% |
| #354 | B*40:01 | C*01:02 | C*07:01 | 8±6 | 100 | 8353 | 5769 | 5500 | nd (02=3) | 2609-7792 | 102-7297 | | |
| #100 | A*11:01 | A*11:01 | B*35:01 | | | | A*24:02 | B*55:01 | C*03:03 | 2 | 27 | 2 | 78% |
| #100 | B*40:02 | C*02:02 | C*04:01 | 5±4 | 100 | 8341 | 8 | 11270 | 14 | 132-169 | 102-11270 | 186-471 | |
| #387 | A*03:01 | A*03:01 | B*07:02 | | | | A*02:05 | B*51:01 | C*16:02 | 15 | 33 | 3 | 87% |
| #387 | B*57:02 | C*07:02 | C*18:01 | 3±3 | 100 | 8180 | nd (01=8498) | 10242 | nd (01=1) | 155-8988 | 132-10242 | 218-568 | |
| #334 | A*01:01 | A*03:01 | B*08:01 | | | | A*23:01 | B*27:02 | C*02:02 | 7 | 28 | 0 | 88% |
| #334 | B*35:03 | C*07:01 | C*07:02 | 4±7 | 100 | 9685 | 8407 | nd (05=5821) | 9 | 105-8978 | 121-8693 | | |
| #188 | A*01:01 | A*03:01 | B*08:01 | | | | A*02:01 | B*51:01 | C*15:02 | 12 | 31 | 1 | 90% |
| #188 | B*37:01 | C*06:02 | C*07:01 | 8±9 | 100 | 11215 | 9506 | 10487 | 47 | 165-11355 | 282-10994 | 135 | |
| #369 | A*26:01 | A*30:01 | B*18:01 | | | | A*02:01 | B*15:01 | C*03:04 | 9 | 32 | 3 | 90% |
| #369 | B*38:01 | C*07:01 | C*12:03 | 4±2 | 100 | 9197 | 2879 | 12245 | 2848 | 172-3194 | 107-12245 | 2848-4268 | |
| #85 | A*02:01 | A*30:01 | B*13:02 | | | | A*68:01 | B*35:01 | C*07:02 | 0 | 11 | 11 | 91% |
| #85 | B*14:02 | C*06:02 | C*08:02 | 8±3 | 100 | 9641 | 18 | 270 | 3383 | | 134-293 | 163-12308 | |
| #177 | A*01:01 | A*02:01 | B*40:02 | | | | A*03:01 | B*07:02 | C*07:02 | 0 | 11 | 10 | 92% |
| #177 | B*57:01 | C*02:02 | C*06:02 | 30±41 | 200 | 9001 | 8 | 3379 | 3045 | | 469-13249 | 606-11695 | |
| #87 | A*01:01 | A*68:01 | B*40:01 | | | | A*02:01 | B*35:01 | C*04:01 | 3 | 28 | 7 | 93% |
| #87 | B*44:02 | C*03:04 | C*07:04 | 2±2 | 100 | 11120 | 358 | 6814 | 183 | 216-358 | 130-6951 | 162-2563 | |
| #376 | A*31:01 | A*32:01 | B*27:05 | | | | A*02:01 | B*51:01 | C*15:02 | 13 | 25 | 0 | 94% |
| #376 | B*40:01 | C*01:02 | C*03:04 | 7±5 | 100 | 8580 | 2947 | 2034 | 42 | 159-3949 | 218-2441 | | |
| #189 | A*24:02 | A*32:01 | B*35:03 | | | | A*01:01 | B*08:01 | C*07:01 | 27 | 3 | 7 | 95% |
| #189 | B*44:29 | C*05:01 | C*12:03 | 1±1 | 100 | 8614 | 11998 | 1800 | nd (02=7107) | 241-13734 | 486-14956 | 256-10851 | |
| #316 | A*02:01 | A*03:01 | B*07:02 | | | | A*24:02 | B*40:01 | C*03:04 | 7 | 39 | 11 | 95% |
| #316 | B*18:01 | C*07:01 | C*07:02 | 67±91 | 400 | 8538 | 8916 | 11888 | 9587 | 3434-11850 | 441-12522 | 539-10443 | |
| #228 | A*01:01 | A*24:02 | B*08:01 | | | | A*02:01 | B*15:01 | C*03:04 | 28 | 42 | 7 | 96% |
| #228 | B*40:02 | C*02:02 | C*07:01 | 117±131 | 500 | 7311 | 10395 | 9061 | 6174 | 814-12334 | 629-13490 | 1120-8945 | |
| #159 | A*01:01 | A*68:01 | B*44:02 | | | | A*02:01 | 6421 | C*06:02 | 4 | 33 | 4 | 96% |
| #159 | B*51:01 | C*14:02 | C*14:03 | 3±3 | 100 | 8709 | 5136 | 6421 | 489 | 165-6456 | 108-6421 | 264-497 | |
| #353 | A*30:02 | A*30:02 | B*18:01 | | | | A*02:01 | B*49:01 | C*07:01 | 12 | 43 | 1 | 97% |
| #353 | B*18:01 | C*05:01 | C*05:01 | 80±133 | 500 | 7792 | 3720 | 9367 | nd (02=1) | 542-9855 | 553-13957 | 824 | |
| #70 | A*11:01 | A*68:01 | B*07:02 | | | | A*02:05 | B*50:01 | C*06:02 | 0 | 46 | 4 | 98% |
| #70 | B*07:02 | C*07:02 | C*07:02 | 15±24 | 100 | 9299 | nd (01=1) | 7891 | 124 | | 135-7891 | 124-288 | |
| #49 | A*02:01 | A*02:01 | B*15:01 | | | | A*01:01 | B*58:01 | C*07:01 | 25 | 15 | 4 | 98% |
| #49 | B*51:01 | C*01:02 | C*03:03 | 33±78 | 300 | 10440 | 11818 | 6305 | nd (02=1418) | 5698-13117 | 455-13108 | 857-6987 | |
| #24 | A*24:02 | A*32:01 | B*07:02 | | | | A*02:01 | B*44:02 | C*05:01 | 26 | 22 | 1 | 99% |
| #24 | B*41:02 | C*07:02 | C*17:03 | 6±9 | 100 | 8627 | 7609 | 6742 | 376 | 529-7928 | 108-12019 | 376 | |
| #275 | A*02:01 | A*02:01 | B*44:05 | | | | A*03:02 | B*18:01 | C*07:01 | 27 | 41 | 8 | 100% |
| #275 | B*57:01 | C*02:02 | C*06:02 | 2±2 | 100 | 8573 | nd (01=8487) | 7916 | nd (02=39) | 157-6798 | 130-7916 | 489-2139 | |
| #384 | A*03:01 | A*33:01 | B*14:02 | | | | A*24:02 | B*40:06 | C*15:02 | 15 | 38 | 12 | 100% |
| #384 | B*15:18 | C*07:04 | C*08:02 | 11±7 | 100 | 8644 | 6122 | 8106 | 8227 | 206-8802 | 107-11983 | 114-8426 | |
| #392 | A*11:01 | A*32:01 | B*35:03 | | | | A*02:01 | B*49:01 | C*07:01 | 10 | 46 | 15 | 100% |
| #392 | B*35:01 | C*04:01 | C*04:01 | 270±96 | 600 | 9291 | 11153 | 7810 | d (02=1257) | 5963-12732 | 772-12611 | 2295-12908 | |

Epitope Analysis Approach

1. Determine antibody specificities to antibody-verified epitopes recorded in the Registry

What proportion reactive alleles carry such epitopes?

2. Analyze the remaining antibody reactivity and look for informative alleles expressing eplet-defined epitopes
3. Such epitopes might be classified as being newly antibody-verified

| Case | Immunizing Allele Antibody-Verified Epitopes on Reactive Alleles | Numbers of Antibody-Verified Epitope-Carrying Reactive Alleles | | | | |
|------|---|---|-------|-------|-------|------|
| | | HLA-A | HLA-B | HLA-C | Total | |
| #186 | 71ATD,161D | 1/1 | 1/7 | 0/2 | 2/10 | 20% |
| #254 | 73TVS,76VS+152RE | 0/0 | 1/2 | 8/8 | 9/10 | 90% |
| #247 | 127K, 62GE | 9/13 | 3/3 | 0/0 | 12/16 | 75% |
| #233 | 41T, 163EW | 1/2 | 17/31 | 0/1 | 18/34 | 53% |
| #373 | 62LQ,66IF+163TEW,76ANT,180E | 7/25 | 8/32 | 0/3 | 15/60 | 25% |
| #129 | 56R,80K,144QL | 3/28 | 2/27 | 7/7 | 12/62 | 19% |
| #355 | 44RT+69TNT,80K | 0/1 | 8/13 | 7/9 | 15/23 | 65% |
| #354 | 56R, 80I,82LR | 8/8 | 19/20 | 0/0 | 27/28 | 96% |
| #100 | 65QIA,69AA,163TEW+65QI | 0/2 | 24/27 | 0/2 | 24/31 | 77% |
| #387 | 44RT+69TNT,127K,144TKH,163LW+65QI | 9/15 | 18/33 | 0/3 | 27/51 | 53% |
| #334 | 65GK,80I,82LR,163EW+73TE | 5/7 | 26/28 | 0/0 | 31/35 | 89% |
| #188 | 80I,131S+163LW,144TKH | 11/12 | 22/31 | 0/1 | 33/44 | 75% |
| #369 | 44RMA,127K,144TKH,163LW,173K | 9/9 | 23/32 | 3/3 | 35/44 | 80% |
| #85 | 44RT,62GE,80K,163TW+65QI,267QE | 0/0 | 1/11 | 2/11 | 3/22 | 14% |
| #177 | 65QIA+76ESN,76VRN,267QE | 0/0 | 11/11 | 10/10 | 21/21 | 100% |
| #87 | 62GE,44RT,80K,131S+163LW | 3/3 | 22/28 | 7/7 | 32/38 | 84% |
| #376 | 131S+163LW,150AAH | 11/13 | 20/25 | 0/0 | 31/38 | 82% |
| #189 | 90D,76ANT,138MI+79GT,267QE | 21/27 | 1/3 | 5/7 | 27/37 | 73% |
| #316 | 41T,73TVS,82LR,166DG,219W | 7/7 | 27/39 | 9/11 | 43/57 | 75% |
| #228 | 44RMA,76VS+152RE,144TKH,144K+76VDT,163LW,173K | 9/28 | 24/42 | 6/7 | 39/77 | 51% |
| #159 | 62GE,76ESN,107W | 4/4 | 31/33 | 0/4 | 35/41 | 85% |
| #353 | 82LR, 41T,131S+163LW,144TKH | 11/12 | 34/43 | 0/1 | 45/56 | 80% |
| #70 | 41T, 69TNT, 131S+163LW | 0/0 | 40/46 | 0/4 | 40/50 | 80% |
| #49 | 65RNA+80I,138MI+79GT,144KR,166DG check c ep | 25/25 | 5/15 | 4/4 | 34/44 | 77% |
| #24 | 62GE,142TKH,80T,166ES | 26/26 | 12/22 | 1/1 | 39/49 | 80% |
| #275 | 73TS, 76ESN, 44RT+69TNT,138MI,144KR,163TEW+65QI | 22/27 | 37/41 | 7/8 | 66/76 | 87% |
| #384 | 21H,41T,80K, 82LR,127K,163EW,167DG | 14/15 | 31/38 | 10/12 | 55/65 | 85% |
| #392 | 1C,41T,62GE,76VRN,82LR+138T,127K,267QE | 10/10 | 27/46 | 15/15 | 52/71 | 73% |

Which additional epitopes might be recognized by antibodies in these first pregnancy sera?

For what proportions of reactive alleles can we identify specific epitopes?

| Case | Additional Epitopes on Reactive Alleles | Total Number of Epitope-Carrying Reactive Alleles | | Numbers of Remaining Reactive Alleles (MFI)* |
|------|--|--|------|---|
| #186 | 6K,71KA, 76ED | 7/10 | 70% | 3 (104-108; C*14:02=691) |
| #254 | None | 9/10 | 90% | 1 (B*73:01=3061) |
| #247 | None | 12/16 | 75% | 4(314-384; A*29:01=650; A*29:02=620) |
| #233 | 80N+66I | 33/34 | 97% | 1 (430) |
| #373 | 62RN+163TW,71STN,245AS | 42/60 | 70% | 18 (117-365;A*03:01;B*37:01=690;B*40:02=831;A*30:01=4947) |
| #129 | 66NV, 80N+45E (TP?)163TEW+182T | 61/62 | 98% | 1 (297) |
| #355 | 21H94I+91G | 18/23 | 78% | 5(243-467;A*33:01=583;B*14:02=1062;B*15:02=1859) |
| #354 | None | 27/28 | 96% | 1 (101) |
| #100 | None | 24/31 | 77% | 7 (102-186;C*03:02=471) |
| #387 | 69TNT+41A,131S+163LW+62R | 41/51 | 80% | 10(132-400; A*33:01=1658;B*54:01=2440) |
| #334 | 71KA | 32/35 | 91% | 3(105-121) |
| #188 | 170RH | 38/44 | 70% | 6 (186-345; B*39:01=447;B*67:01=493; B*15:12=510) |
| #369 | None | 35/44 | 80% | 9 (107-353; B*40:06=602) |
| #85 | 80N+138T,163LW+65QI+70N74Y171Y | 22/22 | 100% | 0 |
| #177 | None | 21/21 | 100% | 0 |
| #87 | 131S | 34/38 | 89% | 4(130-420; B*54:01=750; B*15:12=1732) |
| #376 | 144AH+161E | 33/38 | 87% | 5 (121-343) |
| #189 | 142TKH (TP?) | 33/37 | 86% | 4(156-256; B*15:12=486;C*15:02=606) |
| #316 | 63EI | 49/57 | 86% | 8(441-646;B*15:11=916;B*82:01=1111;B*15:02=1368;B*15:10=1708) |
| #228 | 131S, 166EW+182T | 69/77 | 89% | 8 (297-758;A*80:01=4483) |
| #159 | 80K+14R+16G+138T? | 39/41 | 95% | 2(333; B49:01=5167) |
| #353 | 131S+170RY | 51/56 | 91% | 5(542-718;B*07:02=1310;B*42:01=1639) |
| #70 | 131S | 46/50 | 92% | 4 (124-288) |
| #49 | 138MI,69AA+180Q,166ES (TP) | 43/44 | 98% | 1 (B*15:13=867) |
| #24 | None | 55/65 | 85% | 10 (108-263; B*15:12=8306) |
| #275 | 66NV,69TNT+41A | 76/76 | 100% | 0 |
| #384 | None | 55/65 | 85% | 10 (106-308;B*42:01=404;B*08:01=458;B*82:01=512) |
| #392 | 14R,193PI+183D | 70/71 | 99% | 1 (772) |

*List only alleles with MFI values more than 300 above cut-off values

Conclusion

The HLAMatchmaker version
with antibody-verified class I
epitopes works quite well
but
more antibody-verified epitopes
need to be added


Principles of HLAMatchmaker-based antibody analysis programs

Two examples:

HLA-DP

HLA-C

Luminex Example of Anti-DP Reactivity (Patient Type: DPB1*0201,-)



Anti- HLA-DP Reactivity of Serum Group B Pt 43, Patient types as DPB1*0201,-

| | | | |
|------|---|---|------|
| CON3 | x | x | 76 |
| CON1 | x | x | 107 |
| CON2 | x | x | 121 |
| P 77 | x | x | 7274 |

| | | | | |
|-------|-----------|-----------|------|------|
| P 112 | DPB1*0101 | DPA1*0103 | 2596 | SELF |
| P 113 | DPB1*0101 | DPA1*0201 | 5228 | |
| P 114 | DPB1*0101 | DPA1*0202 | 5932 | |
| P 115 | DPB1*0101 | DPA1*0301 | 2603 | |
| P 116 | DPB1*0201 | DPA1*0103 | 425 | |
| P 117 | DPB1*0301 | DPA1*0103 | 2957 | |
| P 118 | DPB1*0401 | DPA1*0103 | 894 | |
| P 119 | DPB1*0401 | DPA1*0201 | 3807 | |
| P 120 | DPB1*0401 | DPA1*0202 | 3650 | |
| P 121 | DPB1*0401 | DPA1*0301 | 741 | |
| P 122 | DPB1*0402 | DPA1*0103 | 727 | |
| P 123 | DPB1*0402 | DPA1*0301 | 832 | |
| P 124 | DPB1*0501 | DPA1*0103 | 2958 | |
| P 125 | DPB1*0501 | DPA1*0201 | 5863 | |
| P 126 | DPB1*0501 | DPA1*0202 | 4790 | |
| P 127 | DPB1*0501 | DPA1*0301 | 2714 | |
| P 128 | DPB1*1301 | DPA1*0201 | 5080 | |
| P 129 | DPB1*1401 | DPA1*0201 | 5812 | |
| P 130 | DPB1*1701 | DPA1*0201 | 5104 | |
| P 131 | DPB1*1801 | DPA1*0103 | 423 | |
| P 132 | DPB1*1901 | DPA1*0201 | 4878 | |
| P 133 | DPB1*1901 | DPA1*0301 | 2233 | |
| P 134 | DPB1*2801 | DPA1*0202 | 3814 | |

Anti- HLA-DP Reactivity of Serum Group B Pt 43 Patient types as DPB1*0201,-

| | | | |
|------|---|---|------|
| CON3 | x | x | 76 |
| CON1 | x | x | 107 |
| CON2 | x | x | 121 |
| P 77 | x | x | 7274 |

| | | | |
|-------|-----------|-----------|------|
| P 112 | DPB1*0101 | DPA1*0103 | 2596 |
| P 113 | DPB1*0101 | DPA1*0201 | 5228 |
| P 114 | DPB1*0101 | DPA1*0202 | 5932 |
| P 115 | DPB1*0101 | DPA1*0301 | 2603 |
| P 116 | DPB1*0201 | DPA1*0103 | 425 |
| P 117 | DPB1*0301 | DPA1*0103 | 2957 |
| P 118 | DPB1*0401 | DPA1*0103 | 894 |
| P 119 | DPB1*0401 | DPA1*0201 | 3807 |
| P 120 | DPB1*0401 | DPA1*0202 | 3650 |
| P 121 | DPB1*0401 | DPA1*0301 | 741 |
| P 122 | DPB1*0402 | DPA1*0103 | 727 |
| P 123 | DPB1*0402 | DPA1*0301 | 832 |
| P 124 | DPB1*0501 | DPA1*0103 | 2958 |
| P 125 | DPB1*0501 | DPA1*0201 | 5863 |
| P 126 | DPB1*0501 | DPA1*0202 | 4790 |
| P 127 | DPB1*0501 | DPA1*0301 | 2714 |
| P 128 | DPB1*1301 | DPA1*0201 | 5080 |
| P 129 | DPB1*1401 | DPA1*0201 | 5812 |
| P 130 | DPB1*1701 | DPA1*0201 | 5104 |
| P 131 | DPB1*1801 | DPA1*0103 | 423 |
| P 132 | DPB1*1901 | DPA1*0201 | 4878 |
| P 133 | DPB1*1901 | DPA1*0301 | 2233 |
| P 134 | DPB1*2801 | DPA1*0202 | 2814 |

SELF

Questions:

Do the antibodies react with DPB or DPA or both ?

What epitopes are recognized?

Which DP antigens are acceptable mismatches ?

Which DP antigens are unacceptable?

HLAMatchmaker Analysis of Anti- HLA-DP Reactivity of Serum Group B Pt 43

Step 1: Show all eplets in the panel that are mismatched for this patient
(Note the eplets on DPA1*0103 only the “SELF” DPB1*0201)

| | | | | | | | |
|-------|-----------|------------------|------|--------|---|-------------------------------------|--|
| CON3 | x | x | 76 | | | | |
| CON1 | x | x | 107 | | | | |
| CON2 | x | x | 121 | | | | |
| P 77 | x | x | 7274 | mmEp | DPB Eplets | DPA Eplets | |
| P 112 | DPB1*0101 | DPA1*0103 | 2596 | 15 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 113 | DPB1*0101 | DPA1*0201 | 5228 | 15 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 114 | DPB1*0101 | DPA1*0202 | 5932 | 15 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 115 | DPB1*0101 | DPA1*0301 | 2603 | 15 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 116 | DPB1*0201 | DPA1*0103 | 425 | SELF 6 | ,,,,,,,,,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 117 | DPB1*0301 | DPA1*0103 | 2957 | 15 | 8V,11L,,,,,56ED,64DL,65DLK,,70GKR,76V,84DE,87AV,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 118 | DPB1*0401 | DPA1*0103 | 894 | 11 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 119 | DPB1*0401 | DPA1*0201 | 3807 | 11 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 120 | DPB1*0401 | DPA1*0202 | 3650 | 11 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 121 | DPB1*0401 | DPA1*0301 | 741 | 11 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 122 | DPB1*0402 | DPA1*0103 | 727 | 9 | ,,,,,,65DIK,69IKR,70GKR,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 123 | DPB1*0402 | DPA1*0301 | 832 | 9 | ,,,,,,65DIK,69IKR,70GKR,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 124 | DPB1*0501 | DPA1*0103 | 2958 | 13 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 125 | DPB1*0501 | DPA1*0201 | 5863 | 13 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 126 | DPB1*0501 | DPA1*0202 | 4790 | 13 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 127 | DPB1*0501 | DPA1*0301 | 2714 | 13 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 128 | DPB1*1301 | DPA1*0201 | 5080 | 13 | 8V,11L,,,35YA,,56AE,,,,,76I,84DE,87AV,,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 129 | DPB1*1401 | DPA1*0201 | 5812 | 15 | 8V,11L,,,,,56ED,64DL,65DLK,,70GKR,76V,84DE,87AV,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 130 | DPB1*1701 | DPA1*0201 | 5104 | 11 | 8V,11L,,,,,56ED,,,,,,84DE,87AV,,,,,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 131 | DPB1*1801 | DPA1*0103 | 423 | 11 | 8V,,,,,,,,,65DIK,69IKR,70GKR,,84VG,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 132 | DPB1*1901 | DPA1*0201 | 4878 | 10 | ,,,,,,56AE,,,,,76I,84DE,87AV,,,,,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |
| P 133 | DPB1*1901 | DPA1*0301 | 2233 | 10 | ,,,,,,56AE,,,,,76I,84DE,87AV,,,,,,,,,,,,, | 18P,28E,51QA,83T,111K,127L,, | |
| P 134 | DPB1*2801 | DPA1*0202 | 3814 | 11 | ,,,35FA,,,64DL,65DLK,,70GKR,,84VG,,,,,,,,,,,,, | 18P,28E,51RA,83A,111R,127P,, | |

HLAMatchmaker Analysis of Anti- HLA-DP Reactivity of Serum Group B Pt 43

Step 2: After recording DPA1*0103 as negative in the program, identify other negative alleles and enter them in the program

| | | | | | | | |
|-------|-----------|-----------|------|------|------|--|----------------------|
| CON3 | x | x | 76 | | | | |
| CON1 | x | x | 107 | | | | |
| CON2 | x | x | 121 | | | | |
| P 77 | x | x | 7274 | | mmEp | DPB Eplets | DPA Eplets |
| P 112 | DPB1*0101 | DPA1*0103 | 2596 | | 9 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,,, | |
| P 113 | DPB1*0101 | DPA1*0201 | 5228 | | 13 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,,,51RA,83A,111R,127P,, | |
| P 114 | DPB1*0101 | DPA1*0202 | 5932 | | 13 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,,,51RA,83A,111R,127P,, | |
| P 115 | DPB1*0101 | DPA1*0301 | 2603 | | 9 | 8V,,,,35YA,,56AE,,65DIK,69IKR,70GKR,76V,84DE,87AV,,,,,,,,, | |
| P 116 | DPB1*0201 | DPA1*0103 | 425 | SELF | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,, | |
| P 117 | DPB1*0301 | DPA1*0103 | 2957 | | 9 | 8V,11L,,,,,56ED,64DL,65DLK,,70GKR,76V,84DE,87AV,,,,,,,,, | |
| P 118 | DPB1*0401 | DPA1*0103 | 894 | | 5 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,, | |
| P 119 | DPB1*0401 | DPA1*0201 | 3807 | | 9 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 120 | DPB1*0401 | DPA1*0202 | 3650 | | 9 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 121 | DPB1*0401 | DPA1*0301 | 741 | NEG | 5 | ,,,35FA,,56AE,,65DIK,69IKR,70GKR,,,,,,,,, | |
| P 122 | DPB1*0402 | DPA1*0103 | 727 | NEG | 3 | ,,,,,,65DIK,69IKR,70GKR,,,,,,,,, | |
| P 123 | DPB1*0402 | DPA1*0301 | 832 | NEG | 3 | ,,,,,,65DIK,69IKR,70GKR,,,,,,,,, | |
| P 124 | DPB1*0501 | DPA1*0103 | 2958 | | 7 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | |
| P 125 | DPB1*0501 | DPA1*0201 | 5863 | | 11 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 126 | DPB1*0501 | DPA1*0202 | 4790 | | 11 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 127 | DPB1*0501 | DPA1*0301 | 2714 | | 7 | ,,,35LV,,56AE,,65DIK,69IKR,70GKR,,84DE,87AV,,,,,,,,, | |
| P 128 | DPB1*1301 | DPA1*0201 | 5080 | | 11 | 8V,11L,,,35YA,,56AE,,,,,76I,84DE,87AV,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 129 | DPB1*1401 | DPA1*0201 | 5812 | | 13 | 8V,11L,,,,,56ED,64DL,65DLK,,70GKR,76V,84DE,87AV,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 130 | DPB1*1701 | DPA1*0201 | 5104 | | 9 | 8V,11L,,,,,56ED,,,,,84DE,87AV,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 131 | DPB1*1801 | DPA1*0103 | 423 | NEG | 5 | 8V,,,,,,65DIK,69IKR,70GKR,,84VG,,,,,,,,, | |
| P 132 | DPB1*1901 | DPA1*0201 | 4878 | | 8 | ,,,,,,56AE,,,,,76I,84DE,87AV,,,,,,,,, | 51RA,83A,111R,127P,, |
| P 133 | DPB1*1901 | DPA1*0301 | 2233 | | 4 | ,,,,,,56AE,,,,,76I,84DE,87AV,,,,,,,,, | |
| P 134 | DPB1*2801 | DPA1*0202 | 3814 | | 9 | ,,,35FA,,,64DL,65DLK,,70GKR,,84VG,,,,,,,,, | 51RA,83A,111R,127P,, |

HLAMatchmaker Analysis of Anti- HLA-DP Reactivity of Serum Group B Pt 43

Step 3: Remaining eplets after entering the negative DP alleles in the program

| | | | | | | | | |
|-------|-----------|-----------|------|------|---|--------------------------------|--------|--|
| CON3 | x | x | 76 | | | | | |
| CON1 | x | x | 107 | | | | | |
| CON2 | x | x | 121 | | | | | |
| P 77 | x | x | 7274 | mmEp | DPB Eplets | DPA Eplets | | |
| P 112 | DPB1*0101 | DPA1*0103 | 2596 | 4 | ,,,35YA,,,,,,76V,84DE,87AV,,,,,,,,,,,,, | ,,,,,, | | |
| P 113 | DPB1*0101 | DPA1*0201 | 5228 | 8 | ,,,35YA,,,,,,76V,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 114 | DPB1*0101 | DPA1*0202 | 5932 | 8 | ,,,35YA,,,,,,76V,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 115 | DPB1*0101 | DPA1*0301 | 2603 | 4 | ,,,35YA,,,,,,76V,84DE,87AV,,,,,,,,,,,,, | ,,,,,, | | |
| P 116 | DPB1*0201 | DPA1*0103 | 425 | SELF | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,, | |
| P 117 | DPB1*0301 | DPA1*0103 | 2957 | 7 | ,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,,,,,,,,,,,,, | ,,,,,, | | |
| P 118 | DPB1*0401 | DPA1*0103 | 894 | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,, | | |
| P 119 | DPB1*0401 | DPA1*0201 | 3807 | 4 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 120 | DPB1*0401 | DPA1*0202 | 3650 | 4 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 121 | DPB1*0401 | DPA1*0301 | 741 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,, | |
| P 122 | DPB1*0402 | DPA1*0103 | 727 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,, | |
| P 123 | DPB1*0402 | DPA1*0301 | 832 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,, | |
| P 124 | DPB1*0501 | DPA1*0103 | 2958 | 3 | ,,,35LV,,,,,,84DE,87AV,,,,,,,,,,,,, | ,,,,,, | | |
| P 125 | DPB1*0501 | DPA1*0201 | 5863 | 7 | ,,,35LV,,,,,,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 126 | DPB1*0501 | DPA1*0202 | 4790 | 7 | ,,,35LV,,,,,,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 127 | DPB1*0501 | DPA1*0301 | 2714 | 3 | ,,,35LV,,,,,,84DE,87AV,,,,,,,,,,,,, | ,,,,,, | | |
| P 128 | DPB1*1301 | DPA1*0201 | 5080 | 9 | ,11L,,,35YA,,,,,,76I,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 129 | DPB1*1401 | DPA1*0201 | 5812 | 11 | ,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 130 | DPB1*1701 | DPA1*0201 | 5104 | 8 | ,11L,,,,,56ED,,,,,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 131 | DPB1*1801 | DPA1*0103 | 423 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,, | |
| P 132 | DPB1*1901 | DPA1*0201 | 4878 | 7 | ,,,,,,76I,84DE,87AV,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |
| P 133 | DPB1*1901 | DPA1*0301 | 2233 | 3 | ,,,,,,76I,84DE,87AV,,,,,,,,,,,,, | ,,,,,, | | |
| P 134 | DPB1*2801 | DPA1*0202 | 3814 | 6 | ,,,,,,64DL,65DLK,,,,,,,,,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, | | |

HLAMatchmaker Analysis of Anti- HLA-DP Reactivity of Serum Group B Pt 43

Step 4: determine what eplets are shared between reactive alleles

| | | | | | | | |
|-------|-----------|-----------|------|------|------|---|------------------------|
| CON3 | x | x | 76 | | | | |
| CON1 | x | x | 107 | | | | |
| CON2 | x | x | 121 | | | | |
| P 77 | x | x | 7274 | | mmEp | DPB Eplets | DPA Eplets |
| P 112 | DPB1*0101 | DPA1*0103 | 2596 | | 4 | ,,,35YA,,,,,76V,84DE,87AV,,,,, | ,,,, |
| P 113 | DPB1*0101 | DPA1*0201 | 5228 | | 8 | ,,,35YA,,,,,76V,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 114 | DPB1*0101 | DPA1*0202 | 5932 | | 8 | ,,,35YA,,,,,76V,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 115 | DPB1*0101 | DPA1*0301 | 2603 | | 4 | ,,,35YA,,,,,76V,84DE,87AV,,,,, | ,,,, |
| P 116 | DPB1*0201 | DPA1*0103 | 425 | SELF | 0 | ,,,,,,,,,,,,,,,,,,,,,, | ,,,, |
| P 117 | DPB1*0301 | DPA1*0103 | 2957 | | 7 | ,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,,,,, | ,,,, |
| P 118 | DPB1*0401 | DPA1*0103 | 894 | | 0 | ,,,,,,,,,,,,,,,,,,,,,, | ,,,, |
| P 119 | DPB1*0401 | DPA1*0201 | 3807 | | 4 | ,,,,,,,,,,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, |
| P 120 | DPB1*0401 | DPA1*0202 | 3650 | | 4 | ,,,,,,,,,,,,,,,,,,,,,, | ,,51RA,83A,111R,127P,, |
| P 121 | DPB1*0401 | DPA1*0301 | 741 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,, | ,,,, |
| P 122 | DPB1*0402 | DPA1*0103 | 727 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,, | ,,,, |
| P 123 | DPB1*0402 | DPA1*0301 | 832 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,, | ,,,, |
| P 124 | DPB1*0501 | DPA1*0103 | 2958 | | 3 | ,,,35LV,,,,,84DE,87AV,,,,, | ,,,, |
| P 125 | DPB1*0501 | DPA1*0201 | 5863 | | 7 | ,,,35LV,,,,,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 126 | DPB1*0501 | DPA1*0202 | 4790 | | 7 | ,,,35LV,,,,,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 127 | DPB1*0501 | DPA1*0301 | 2714 | | 3 | ,,,35LV,,,,,84DE,87AV,,,,, | ,,,, |
| P 128 | DPB1*1301 | DPA1*0201 | 5080 | | 9 | ,11L,,,35YA,,,,,76I,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 129 | DPB1*1401 | DPA1*0201 | 5812 | | 11 | ,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 130 | DPB1*1701 | DPA1*0201 | 5104 | | 8 | ,11L,,,,,56ED,,,,,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 131 | DPB1*1801 | DPA1*0103 | 423 | NEG | 0 | ,,,,,,,,,,,,,,,,,,,,,, | ,,,, |
| P 132 | DPB1*1901 | DPA1*0201 | 4878 | | 7 | ,,,,,,,,,,76I,84DE,87AV,,,,, | ,,51RA,83A,111R,127P,, |
| P 133 | DPB1*1901 | DPA1*0301 | 2233 | | 3 | ,,,,,,,,,,76I,84DE,87AV,,,,, | ,,,, |
| P 134 | DPB1*2801 | DPA1*0202 | 3814 | | 6 | ,,,,,,64DL,65DLK,,,,, | ,,51RA,83A,111R,127P,, |


HLAMatchmaker-Predicted Unacceptable DPB Alleles not Present in the Luminex Panel

| | | | |
|-----------|---|-----------|--|
| DPB1*0601 | 8V,11L,,,,,56ED,64DL,65DLE,,,,84DE,87AV,, | DPB1*5601 | 8V,11L,,,,,64DL,65DLK,,,76V,84DE,87AV,, |
| DPB1*0801 | ,,,,,,76V,84DE,87AV,, | DPB1*5701 | ,,,,,56ED,64DL,65DLK,,,76V,84DE,,, |
| DPB1*1601 | ,,,,,,84DE,87AV,, | DPB1*5801 | 8V,11L,,,35LV,,,,,,84DE,87AV,, |
| DPB1*2001 | 8V,11L,,,,,56ED,64DL,65DLK,,,84DE,87AV,, | DPB1*6301 | ,,,35LV,,,,,,84DE,87AV,, |
| DPB1*2101 | 8V,11L,,,35LV,,,,,,84DE,87AV,, | DPB1*6501 | ,,,35YA,,,,,,76V,84DE,87AV,, |
| DPB1*2201 | ,,,35LV,,,,,,84DE,87AV,, | DPB1*6701 | 8V,11L,,,,,64DL,65DLK,,,76V,84DE,87AV,, |
| DPB1*2501 | 8V,11L,,,,,64DL,65DLK,,,76V,84DE,87AV,, | DPB1*6801 | ,,,,,,76V,84DE,87AV,, |
| DPB1*2601 | 8V,11L,,,35YA,,,,,,76V,84DE,87AV,, | DPB1*6901 | 8V,11L,,,,,56ED,64DL,65DLR,,,84DE,87AV,, |
| DPB1*2701 | 8V,11L,,,35YA,,,,,,84DE,87AV,, | DPB1*7001 | 8V,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,, |
| DPB1*2901 | 8V,11L,,,,,56ED,64DL,65DLE,,,76V,84DE,87AV,, | DPB1*7601 | 8V,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,, |
| DPB1*3001 | 8V,11L,,,,,,84DE,87AV,, | DPB1*7801 | 8V,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,, |
| DPB1*3101 | ,,,,,64DL,65DLK,,,84DE,87AV,, | DPB1*7901 | 8V,11L,,,,,,76V,84DE,87AV,, |
| DPB1*3501 | 8V,11L,,,,,56ED,,,,,76V,84DE,87AV,, | DPB1*8401 | ,,,,,,76V,84DE,87AV,, |
| DPB1*3601 | 8V,11L,,,35LV,,,,,,84DE,87AV,, | DPB1*8501 | 8V,11L,,,35YA,,,,,,84DE,87AV,91H, |
| DPB1*3701 | 8V,11L,,,,,,76V,84DE,87AV,, | DPB1*8701 | 8V,11L,,,,,64DL,65DLK,,,84DE,87AV,, |
| DPB1*3801 | ,,,35LV,,,,,,84DE,87AV,, | DPB1*8801 | 8V,11L,,,,,56ED,,,,,76V,84DE,87AV,, |
| DPB1*4401 | 8V,11L,,,35LV,56ED,64DL,65DLE,,,76V,84DE,87AV,, | DPB1*8901 | 8V,,,,,35YA,,,,,,84DE,87AV,, |
| DPB1*4501 | 8V,11L,,,,,64DL,65DLK,,,76V,84DE,87AV,, | DPB1*9001 | 8V,,,,,,76V,84DE,87AV,, |
| DPB1*5001 | 8V,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,, | DPB1*9101 | 8V,11L,,,,,56ED,64DL,65DLK,,,84DE,87AV,, |
| DPB1*5201 | 8V,11L,,,,,64DL,65DLK,,,76V,84DE,87AV,, | DPB1*9201 | 8V,11L,,,,,56ED,64DL,65DLK,,,76V,84DE,87AV,, |
| DPB1*5401 | 8V,11L,,,,,,76V,84DE,87AV,, | DPB1*9301 | 8V,11L,,,,,,84DE,87AV,, |
| DPB1*5501 | 8V,11L,,,,,,84DE,87AV,, | DPB1*9701 | ,,,35LV,43WR,,,,,,84DE,87AV,, |
| | | DPB1*9801 | 8V,11L,,,,,56ED,,,,,84DE,87AV,, |

HLAMatchmaker-Predicted Acceptable DPB Alleles not Present in the Luminex Panel

| | | | | | |
|-----------|---|--|-----------|---|---|
| DPB1*0202 | 1 | ,,,35LV,,,,,,,,, | DPB1*6001 | 1 | ,,,,,64DN,,,,,,,,, |
| DPB1*1501 | 6 | 8V,,,33Q,35YA,,,64DL,65DLR,,,,,84VG,,, | DPB1*6201 | 2 | ,,,35LV,,,,,,,,,84VG,,, |
| DPB1*1801 | 2 | 8V,,,,,,,,,,84VG,,, | DPB1*6601 | 2 | 8V,11L,,,,,,,,,,,,, |
| DPB1*2301 | 0 | ,,,,,,,,,,,,, | DPB1*7101 | 0 | ,,,,,,,,,,,,, |
| DPB1*2401 | 0 | ,,,,,,,,,,,,, | DPB1*7201 | 2 | ,,,,,64DL,65DLK,,,,,,,, |
| DPB1*2801 | 3 | ,,,,,,64DL,65DLK,,,84VG,,, | DPB1*7301 | 3 | ,,,,,,64DL,65DLK,,,76V,,, |
| DPB1*3201 | 1 | ,,,,,56EV,,,,,,,,, | DPB1*7401 | 7 | 8V,11L,,33Q,35YA,,,64DL,65DLR,,,,,84VG,,, |
| DPB1*3301 | 0 | ,,,,,,,,,,,,, | DPB1*7501 | 1 | ,,,,,,76V,,,,, |
| DPB1*3401 | 4 | ,,,35LV,,,64DL,65DLK,,,84VG,,, | DPB1*7701 | 0 | ,,,,,,,,,,,,, |
| DPB1*3901 | 1 | ,,,35YA,,,,,,,,,, | DPB1*8001 | 1 | ,,,,,56ED,,,,,,,,, |
| DPB1*4001 | 2 | ,,,35YA,,,,,,,,,84VG,,, | DPB1*8101 | 0 | ,,,,,,,,,,,,, |
| DPB1*4101 | 3 | ,,,,,,64DF,65DFE,69FER,,,,,,,, | DPB1*8201 | 1 | ,,,,,,,,,,,,,91H, |
| DPB1*4601 | 1 | ,,,,,56ED,,,,,,,,, | DPB1*8301 | 2 | ,,,,,,64DF,65DFK,,,,,,,, |
| DPB1*4701 | 0 | ,,,,,,,,,,,,, | DPB1*8601 | 3 | 8V,11L,,,,,56ED,,,,,,,,, |
| DPB1*4801 | 1 | ,,,35LV,,,,,,,,, | DPB1*9401 | 1 | ,,,43GW,,,,,,,,, |
| DPB1*4901 | 1 | ,,,35YA,,,,,,,,, | DPB1*9501 | 3 | ,,,35LV,,,64DL,65DLE,,,,,,,, |
| DPB1*5101 | 0 | ,,,,,,,,,,,,, | DPB1*9601 | 2 | ,,,35YA,,,64HI,,,,,,,,, |
| DPB1*5301 | 2 | ,,,35YA,,,,,,,,,84VG,,, | DPB1*9901 | 1 | ,,28D,,,,,,,,,,,,, |

Molecular Locations of Positions 84 and 87



84DE+87AV Corresponds to an Antibody-Defined Epitope

HLA-DP Epitope Typing Using Monoclonal Antibodies

William H. Marshall, Sheila Drover, Dianne Codner, Jane Gamberg, M. Douglas Copp, Hong-Wei Liu, Lang-Tuo Deng, and H. Banfield Younghusband

Human Immunology 59, 189–197 (1998)

TABLE 4 Proposed specificities of the monoclonal antibodies

| | | |
|-----------|-------|-----------------|
| NFLD.M66 | IgG1 | DPB1:84–87:DEAV |
| NFLD.M126 | IgG1 | DPB1:84–87:DEAV |
| NFLD.M120 | IgG2a | DPB1:85–87:GPM |
| NFLD.M123 | IgM | DPB1:84–87:GGPM |

Frequency of Anti-84DE,87AV Reactive Antibodies

Our experience:

- 78/355 (22%) HLA class II-reactive patient sera have anti-HLA-DP antibodies
- 33/78 (42%) react with 84DE,87AV

Frequency of Anti-84DE,87AV Reactive Antibodies

Our experience:

- 78/355 (22%) HLA-reactive patient sera have anti-HLA-DP antibodies
- 33/78 (42%) react with 84DE,87AV

Danny Youngs (Seattle) ASHI Quarterly 2004

29/67 (43%) of patients with anti-DP antibodies
react with DEAV at position 84-87

HLAMatchmaker Determination of Unacceptable DPA Alleles

| | Eplet | Unacceptable Alleles |
|-----|-------|---------------------------|
| DPA | 51RA | *0108*0201*0202*0203*0401 |
| DPA | 83A | *0105*0201*0202*0203*0401 |
| DPA | 111R | *0201*0202 |
| DPA | 127P | *0201*0202*0401 |

HLA-C Antibodies

Example of Luminex Screen for Antibodies Specific for HLA-C (1)

| | |
|-----------|------------|
| Patient | HLA-C Type |
| Patient | Cw*0602 |
| Patient | Cw*0702 |
| Imm Donor | Cw*0102 |
| Imm Donor | Cw*0701 |

| Bead | Allele | MFI | RX |
|------|---------|-------|------|
| 1 | NegCont | 880 | |
| 2 | PosCont | 12591 | |
| 83 | Cw*0102 | 9850 | IMM |
| 84 | Cw*0202 | 1054 | |
| 98 | Cw*0302 | 2118 | |
| 97 | Cw*0303 | 5498 | |
| 99 | Cw*0304 | 4991 | |
| 85 | Cw*0401 | 1393 | |
| 86 | Cw*0501 | 1330 | |
| 87 | Cw*0602 | 813 | SELF |
| 88 | Cw*0702 | 624 | SELF |
| 89 | Cw*0801 | 2465 | |
| 90 | Cw*1203 | 767 | |
| 91 | Cw*1402 | 7769 | |
| 92 | Cw*1502 | 738 | |
| 93 | Cw*1601 | 4305 | |
| 94 | Cw*1701 | 798 | |
| 95 | Cw*1802 | 1099 | |

What is the cut-off point to distinguish between positive and negative reactions?

Example of Luminex Screen for Antibodies Specific for HLA-C (2)

| | |
|-----------|------------|
| Patient | HLA-C Type |
| Patient | Cw*0602 |
| Patient | Cw*0702 |
| Imm Donor | Cw*0102 |
| Imm Donor | Cw*0701 |

| Cut Off | | 1400 | |
|---------|---------|-------|------|
| Bead | Allele | MFI | RX |
| 1 | NegCont | 880 | |
| 2 | PosCont | 12591 | |
| 83 | Cw*0102 | 9850 | IMM |
| 84 | Cw*0202 | 1054 | NEG |
| 98 | Cw*0302 | 2118 | POS |
| 97 | Cw*0303 | 5498 | POS |
| 99 | Cw*0304 | 4991 | POS |
| 85 | Cw*0401 | 1393 | NEG |
| 86 | Cw*0501 | 1330 | NEG |
| 87 | Cw*0602 | 813 | SELF |
| 88 | Cw*0702 | 624 | SELF |
| 89 | Cw*0801 | 2465 | POS |
| 90 | Cw*1203 | 767 | NEG |
| 91 | Cw*1402 | 7769 | POS |
| 92 | Cw*1502 | 738 | NEG |
| 93 | Cw*1601 | 4305 | POS |
| 94 | Cw*1701 | 798 | NEG |
| 95 | Cw*1802 | 1099 | NEG |

Antibody react with Cw*0102 of donor
 Reactivity with donor's Cw*0701? No allele in panel
 What HLA-C epitopes are recognized?

Do HLAMatchmaker analysis
 What are the mismatched eplets?

Example of Luminex Screen for Antibodies Specific for HLA-C (3)

| | |
|-----------|--|
| | HLA-C Type |
| Patient | Cw*0602 |
| Patient | Cw*0702 |
| Imm Donor | Cw*0102 ,1C,6K,,,,,,69QRQT,,,,,77TVS,,,,,,,,,,,,,219W,,,,, |
| Imm Donor | Cw*0701 ,,,,,,,,,,65QNR,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |

| | Cut Off | 1400 | | |
|------|---------|-------|------|--|
| Bead | Allele | MFI | RX | Mismatched Eplets |
| 1 | NegCont | 880 | | |
| 2 | PosCont | 12591 | | |
| 83 | Cw*0102 | 9850 | IMM | ,1C,6K,,,,,69QRQT,,,,,77TVS,,,,,219W,,,,, |
| 84 | Cw*0202 | 1054 | NEG | ,,21H,,,,,69QRQT,,,,,77TVN,,,,,211T,,,,, |
| 98 | Cw*0302 | 2118 | POS | ,,21H,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,163LW,,173K,,,219W,,,,, |
| 97 | Cw*0303 | 5498 | POS | ,,21H,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,163LW,,173K,,,219W,,,,, |
| 99 | Cw*0304 | 4991 | POS | ,,21H,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,163LW,,173K,,,219W,,,,, |
| 85 | Cw*0401 | 1393 | NEG | ,,17WR,,,,,113YN,116F,,,,,219W,,,275KP,, |
| 86 | Cw*0501 | 1330 | NEG | ,,35Q,,,46QGE,,,,,69QRQT,,,,,77TVN,,,,,113YN,116F,138K,,,,,177KT,,,,,275GP,, |
| 87 | Cw*0602 | 813 | SELF | |
| 88 | Cw*0702 | 624 | SELF | |
| 89 | Cw*0801 | 2465 | POS | ,,35Q,,,46QGE,,,,,69QRQT,,,,,77TVS,,,,,113YN,116F,,,,,152RT,,177KT,,,,,275GP,, |
| 90 | Cw*1203 | 767 | NEG | |
| 91 | Cw*1402 | 7769 | POS | ,,,,,69QRQT,,,,,77TVS,,,,,219W,,,,, |
| 92 | Cw*1502 | 738 | NEG | ,,21H,,,,,65QNR,,69QRQT,,,,,77TVN,,,,,94II,,,,, |
| 93 | Cw*1601 | 4305 | POS | ,,,,,69QRQT,,,,,77TVS,,,,,193LV,,,,, |
| 94 | Cw*1701 | 798 | NEG | ,,,,,113YN,116F,143ISQ,,,,,184R,,,270C,,,275KP,, |
| 95 | Cw*1802 | 1099 | NEG | ,,,,,113YN,116F,,,,,219W,,,275KP,, |

These are the eplets on HLA-C alleles of donor and the panel

Cw*1203 has no mismatched eplets for this patient

One or more eplets on donor's Cw*0102 must react with antibody

Any reactivity with 65QNR of donor's Cw*0701?

Example of Luminex Screen for Antibodies Specific for HLA-C (4)

| | HLA-C Type |
|-----------|---|
| Patient | Cw*0602 |
| Patient | Cw*0702 |
| Imm Donor | Cw*0102 ,1C,6K,,,,,,69QRQT,,,,77TVS,,,,,,,,,,,,,219W,,,,, |
| Imm Donor | Cw*0701 ,,,,,,65QNR,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |

| | Cut Off | 1400 | | | |
|------|---------|-------|------|--|--|
| Bead | Allele | MFI | RX | Mismatched Eplets | |
| 1 | NegCont | 880 | | ,, | |
| 2 | PosCont | 12591 | | ,, | |
| 83 | Cw*0102 | 9850 | IMM | ,1C,6K,,,,,,69QRQT,,,,,77TVS,,,,,,,,,,,,,219W,,,,, | |
| 84 | Cw*0202 | 1054 | NEG | ,21H,,,,,,69QRQT,,,,,77TVN,,,,,,,,,,,,,211T,,,,, | |
| 98 | Cw*0302 | 2118 | POS | ,21H,,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,,,,,,163LW,173K,,219W,,,,, | |
| 97 | Cw*0303 | 5498 | POS | ,21H,,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,,,,,,163LW,173K,,219W,,,,, | |
| 99 | Cw*0304 | 4991 | POS | ,21H,,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,,,,,,163LW,173K,,219W,,,,, | |
| 85 | Cw*0401 | 1393 | NEG | ,17WR,,,,,,,,,,,,,113YN,116F,,,,,,,,,,219W,,275KP,, | |
| 86 | Cw*0501 | 1330 | NEG | ,35Q,,,46QGE,,,69QRQT,,,,,77TVN,,,,,,,,,,113YN,116F,138K,,,,,177KT,,,,,275GP,, | |
| 87 | Cw*0602 | 813 | SELF | ,, | |
| 88 | Cw*0702 | 624 | SELF | ,, | |
| 89 | Cw*0801 | 2465 | POS | ,35Q,,,46QGE,,,69QRQT,,,,,77TVS,,,,,,,,,,113YN,116F,,,,,152RT,,177KT,,,,,275GP,, | |
| 90 | Cw*1203 | 767 | NEG | ,, | |
| 91 | Cw*1402 | 7769 | POS | ,,,,,,69QRQT,,,,,77TVS,,,,,,,,,,,,,219W,,,,, | |
| 92 | Cw*1502 | 738 | NEG | ,21H,,,,, <u>65QNR</u> ,69QRQT,,,,,77TVN,,,,,94II,,,,,,,,,,,,, | |
| 93 | Cw*1601 | 4305 | POS | ,,,,,,69QRQT,,,,,77TVS,,,,,,,,,,,,,193LV,,,,, | |
| 94 | Cw*1701 | 798 | NEG | ,,,,,,,,,,,,,113YN,116F,143ISQ,,,,,,,,,,184R,,270C,,275KP,, | |
| 95 | Cw*1802 | 1099 | NEG | ,,,,,,,,,,,,,113YN,116F,,,,,,,,,,219W,,275KP,, | |

No reactivity with 65QNR of donor's Cw*0701

No reactivity with the eplets of the negative alleles

Cw*0202, Cw*0401, Cw*0501, Cw*1502, Cw*1701 and Cw*1801

Record these negative alleles in the HLA Matchmaker program

| | HLA-C Type |
|-----------|--|
| Patient | Cw*0602 |
| Patient | Cw*0702 |
| Imm Donor | Cw*0102,1C,6K,,,,,69QRQT,,,77TVS,,,,,,,,,,,,,,,,,,,,,219W,,,,, |
| Imm Donor | Cw*0701,,,,,,,,,65QNR,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |

| | Cut Off | 1400 | | | | | |
|------|---------|-------|------|--|--|--|--|
| Bead | Allele | MFI | RX | Mismatched Eplets | | Eplets on Reactive Alleles | |
| 1 | NegCont | 880 | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | |
| 2 | PosCont | 12591 | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | |
| 83 | Cw*0102 | 9850 | IMM | ,1C,6K,,,,,69QRQT,,77TVS,,,,,,,,,,,,,219W,,, | | ,1C,6K,,,,,,,,,77TVS,,,,,,,,,,,,, | |
| 84 | Cw*0202 | 1054 | NEG | ,21H,,,,,,,,,69QRQT,,77TVN,,,,,,,,,,,,,211T,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 98 | Cw*0302 | 2118 | POS | ,21H,,,,,,,,,69QRQT,,77TVS,,,,,94II,,,,,,,,,163LW,173K,,,219W,,, | | ,,,,,,,,,77TVS,,,,,,,,,,,,,163LW,173K,,, | |
| 97 | Cw*0303 | 5498 | POS | ,21H,,,,,,,,,69QRQT,,77TVS,,,,,94II,,,,,,,,,163LW,173K,,,219W,,, | | ,,,,,,,,,77TVS,,,,,,,,,,,,,163LW,173K,,, | |
| 99 | Cw*0304 | 4991 | POS | ,21H,,,,,,,,,69QRQT,,77TVS,,,,,94II,,,,,,,,,163LW,173K,,,219W,,, | | ,,,,,,,,,77TVS,,,,,,,,,,,,,163LW,173K,,, | |
| 85 | Cw*0401 | 1393 | NEG | ,17WR,,,,,,,,,,,,,,,,,113YN,116F,,,,,,,,,219W,,,275KP,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 86 | Cw*0501 | 1330 | NEG | ,35Q,,,46QGE,,,69QRQT,,77TVN,,,,,,,,,113YN,116F,138K,,,,,177KT,,,,,2i | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 87 | Cw*0602 | 813 | SELF | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 88 | Cw*0702 | 624 | SELF | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 89 | Cw*0801 | 2465 | POS | ,35Q,,,46QGE,,,69QRQT,,77TVS,,,,,,,,,113YN,116F,,,,,152RT,,177KT,,,,,; | | ,,,,,,,,,77TVS,,,,,,,,,,,,,152RT | |
| 90 | Cw*1203 | 767 | NEG | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 91 | Cw*1402 | 7769 | POS | ,,,,,,,,,69QRQT,,77TVS,,,,,,,,,,,,,219W,,, | | ,,,,,,,,,77TVS,,,,,,,,,,,,, | |
| 92 | Cw*1502 | 738 | NEG | ,21H,,,,,,,,,65QNR,,69QRQT,,77TVN,,,,,94II,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 93 | Cw*1601 | 4305 | POS | ,,,,,,,,,69QRQT,,77TVS,,,,,,,,,,,,,193LV,,, | | ,,,,,,,,,77TVS,,,,,,,,,,,,,193LV,,, | |
| 94 | Cw*1701 | 798 | NEG | ,,,,,,,,,,,,,,,,,113YN,116F,143ISQ,,,,,,184R,,270C,,,275KP,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| 95 | Cw*1802 | 1099 | NEG | ,,,,,,,,,,,,,,,,,113YN,116F,,,,,,,,,219W,,,275KP,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |

The 77TVS eplet of donor's Cw*0102 is present on all reactive alleles
77TVS is an unacceptable mismatch

Example of Luminex Screen for Antibodies Specific for HLA-C (6)

| | HLA-C Type |
|---------|------------|
| Patient | Cw*0602 |
| Patient | Cw*0702 |

Imm Donor Cw*0102,1C,6K,,,,,69QRQT,,,,**77TVS**,,,,,,,,,,,,,219W,,,,,
Imm Donor Cw*0701,,,,,65QNR,,,,,,,,,,,,,,,,,,,,,,,,,,,,,

| Cut Off | | 1400 | | | | |
|---------|---------|-------|------|--|--|----------------------------|
| Bead | Allele | MFI | RX | Mismatched Eplets | | Eplets on Reactive Alleles |
| 1 | NegCont | 880 | | | | |
| 2 | PosCont | 12591 | | | | |
| 83 | Cw*0102 | 9850 | IMM | ,1C,6K,,,,,69QRQT,,,,,77TVS,,,,,219W,,,,, | | ,1C,6K,,,,,77TVS,,,,, |
| 84 | Cw*0202 | 1054 | NEG | ,,21H,,,,,69QRQT,,,,,77TVN,,,,,211T,,,,, | | |
| 98 | Cw*0302 | 2118 | POS | ,,21H,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,163LW,173K,,219W,,,,, | | 77TVS,,,,,163LW,173K,, |
| 97 | Cw*0303 | 5498 | POS | ,,21H,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,163LW,173K,,219W,,,,, | | 77TVS,,,,,163LW,173K,, |
| 99 | Cw*0304 | 4991 | POS | ,,21H,,,,,69QRQT,,,,,77TVS,,,,,94II,,,,,163LW,173K,,219W,,,,, | | 77TVS,,,,,163LW,173K,, |
| 85 | Cw*0401 | 1393 | NEG | ,,17WR,,,,,113YN,116F,,,,,219W,,,275KP,, | | |
| 86 | Cw*0501 | 1330 | NEG | ,,35Q,,,46QGE,,,,,69QRQT,,,,,77TVN,,,,,113YN,116F,138K,,,,,177KT,,,,,219W,,,,, | | |
| 87 | Cw*0602 | 813 | SELF | | | |
| 88 | Cw*0702 | 624 | SELF | | | |
| 89 | Cw*0801 | 2465 | POS | ,,35Q,,,46QGE,,,,,69QRQT,,,,,77TVS,,,,,113YN,116F,,,,,152RT,,177KT,,,,, | | 77TVS,,,,,152RT,,,,, |
| 90 | Cw*1203 | 767 | NEG | | | |
| 91 | Cw*1402 | 7769 | POS | ,,,,,69QRQT,,,,,77TVS,,,,,219W,,,,, | | 77TVS,,,,, |
| 92 | Cw*1502 | 738 | NEG | ,,21H,,,,,65QNR,,69QRQT,,,,,77TVN,,,,,94II,,,,, | | |
| 93 | Cw*1601 | 4305 | POS | ,,,,,69QRQT,,,,,77TVS,,,,,193LV,,,,, | | 77TVS,,,,,193LV,,,,, |
| 94 | Cw*1701 | 798 | NEG | ,,,,,113YN,116F,143ISQ,,,,,184R,,270C,,,275KP,, | | |
| 95 | Cw*1802 | 1099 | NEG | ,,,,,113YN,116F,,,,,219W,,,275KP,, | | |

Acceptable Alleles

Cw*0202*0401*0403*0404*0405*0501*0602*0603*0604*0701*0702*0704*0706*0707*0710*1202*1203*1204*1502*1503*1505*1506*1509*1701*1801

Eplet Unacceptable Alleles

77TVS B*4601Cw*0102*0103*0302*0303*0304*0305*0306*0310*0801*0802*0803*0804*1402*1403*1601*1604