

HLA Epitopes-Special Seminar

Sydney, Australia

March 5-6, 2015

HLA Antibodies in Transplantation

- HLA antibodies cause allograft rejection and transplant failure
- HLA antibodies recognize epitopes

Therefore

HLA epitopes are important in transplantation

Antibody Responses to HLA Affect Transplant Outcome

- Class I HLA antigens
 - HLA-A and HLA-B
 - HLA-C
 - MICA
- Class II HLA antigens
 - DRB1
 - DRB3, 4, 5
 - DQB and DQA
 - DPB and DPA

HLA Antigens Have Multiple Epitopes

- “Private” epitopes and “public” epitopes shared between cross-reacting antigens (CREGs)
- Epitopes defined by polymorphic amino acid residues
- Structurally based matching with HLAMatchmaker
 - Original version (2002): “triplets”, i.e. linear three-residue sequences
 - New version (2006): “eplets”, i.e. patches of residues within a 3 Angstrom radius of a polymorphic residue on the molecular surface.

What Do We Need to Know about HLA Epitopes?

- Structure and topography of epitopes on mismatched HLA antigens
- Structural interactions of HLA antigen-antibody complexes
- Epitope antigenicity: reactivity with antibody *in vivo* and in different assays
- Epitope immunogenicity: ability to induce specific antibodies

Objectives

- Structural concepts of HLA epitopes
 - Triplets, eplets, HLAMatchmaker, structural epitopes
 - Website-based International Registry of HLA-A,B,C, DR, DQ, DP and MICA epitopes
- HLA epitope antigenicity
 - HLAMatchmaker analysis of sera with antibodies
 - Epitope specificities of HLA antibodies tested in different assays
 - Mismatch acceptability for sensitized patients

Objectives (continued)

- HLA epitope immunogenicity
 - Epitope loads and HLA antibody responses
 - Physiochemical characteristics of HLA epitopes
 - Nonself-self paradigm
 - Mismatch permissibility for non-sensitized patients
- Clinical relevance of HLA matching at the epitope level
 - Solid organ transplantation
 - Platelet transfusion
 - Stem cell transplantation
- Future plans for HLA epitopes

Duquesnoy's Presentations

- More than 350 Powerpoint slides
- References to almost 100 publications about HLA epitopes since 2002 (List will be provided)

Lecture 1

Basic Aspects of HLA Epitope Structure

How Do Antibodies React with HLA Molecules?

Consider the structure of

- HLA molecules
- Antibodies
- Epitope-Paratope Interphase



What do you see?


Fotoğraf: George Steinmetz

© 2005 National Geographic Society. Her hakkı saklıdır.


Dev Develer

National Geographic Türkiye, Şubat 2005

Binding of Antibody to an Epitope




Antigen-Antibody Complex




How Do Antibodies Bind to Epitopes on Proteins?

CDRs of Antibody



- Complementarity Determining Regions H1, H2, H3, on heavy chain and L1, L2 and L3 on light chain represent the “binding face” or “paratope” of an antibody
- The antibody specificity is often determined by a centrally located loop (CDR-H3) that binds to a short amino acid sequence of the protein antigen
- Other CDRs serve as contact sites to stabilize binding to antigen (they play a role in the affinity of antibody)

Antigen-Antibody Complex



MMDB Database

<http://www.ncbi.nih.gov/Structure>

- Large database of crystal structures
 - Complexes of protein antigens and antibody domains
 - Class I and class II HLA molecules
 - MHC-T-Cell Receptor complexes
- **Cn3D molecular viewer**
 - “Space-fill” command identifies surface-exposed amino acid residues
 - “Select by distance” command permits an assessment of epitope size and the distance between residues
- Basic Local Alignment Search Tool (BLAST)
 - Identifies sequence differences between antigen proteins and homologous self-proteins of the antibody producer

Commonly Studied Protein Antigens

- Hen Egg Lysozyme (HEL)
- Horse Cytochrome C
- Human Interferon- γ Receptor

With HLAMatchmaker on my mind I have analyzed data in more than 100 publications

Mouse Monoclonal Antibodies Against Hen Egg Lysozyme

What are the differences between
mouse and hen egg lysozyme?

Sequence Differences between HEL and Mouse Lysozymes

intestinal milk HEL	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	
	K	V	Y	N	R	C	E	L	A	R	I	L	K	R	N	G	M	D	G	Y	R	G	V	K	L	
	K	V	Y	E	R	C	E	F	A	R	T	L	K	R	N	G	M	A	G	Y	Y	G	V	S	L	
intestinal milk HEL	-	-	F	G	-	-	-	-	-	A	A	M	-	-	H	-	L	-	N	-	-	-	Y	S	-	
	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	
	A	D	W	V	C	L	A	Q	H	E	S	N	Y	N	T	R	A	T	N	Y	N	R	D	R	S	
intestinal milk HEL	A	D	W	V	C	L	A	Q	H	E	S	N	Y	N	T	R	A	T	N	Y	N	R	D	Q	S	
	G	N	-	-	-	A	-	K	F	-	-	-	F	-	-	Q	-	-	-	R	-	T	-	G	-	
	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	
intestinal milk HEL	T	D	Y	G	I	F	Q	I	N	S	R	Y	W	C	N	D	G	K	T	P	R	S	K	N	A	
	T	D	Y	G	I	F	Q	I	N	S	R	Y	W	C	N	D	G	K	T	P	R	A	V	N	A	
	-	-	-	-	-	L	-	-	-	-	-	W	-	-	-	-	-	R	-	-	G	-	R	-	L	
intestinal milk HEL	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100	
	C	G	I	N	C	S	A	L	L	Q	D	D	I	T	A	A	I	Q	C	A	K	R	V	V	R	
	C	G	I	N	C	S	A	L	L	Q	D	D	I	T	A	A	I	Q	C	A	K	R	V	V	R	
intestinal milk HEL	-	N	-	P	-	-	-	-	-	S	S	-	-	-	-	S	V	N	-	-	-	K	I	-	S	
	101	102	103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126
	D	P	Q	G	I	R	A	W	V	A	W	R	T	Q	C	Q	N	R	D	L	S	Q	Y	I	R	N
intestinal milk HEL	D	P	Q	G	I	R	A	W	V	A	W	R	A	H	C	Q	N	R	D	L	S	Q	Y	I	R	N
	-	G	N	-	M	N	-	-	-	-	-	-	N	R	-	K	G	T	-	V	Q	A	W	-	-	G

Which Residues are Critical for Antigen-Antibody Contact?

- Antibody reactivity with natural variants

Genetic Marker Panel 1 (28 Markers)																											
Sample ID	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27
Mouse intest	K	V	Y	N	R	C	E	L	A	R	I	L	K	R	N	G	M	D	G	Y	R	G	V	K	L	A	D
Mouse milk	K	V	Y	E	R	C	E	F	A	R	T	L	K	R	N	G	M	A	G	Y	Y	G	V	S	L	A	D
HEL	-	-	F	G	-	-	-	-	-	A	A	M	-	-	H	-	L	-	N	-	-	-	Y	S	-	G	N
Cal Quail	-	-	F	G	-	-	-	-	-	A	A	M	-	-	H	-	L	-	N	-	-	-	Y	S	-	G	N
BW Quail	-	-	F	G	-	-	-	-	-	A	A	M	-	-	H	-	L	-	N	-	-	-	Y	S	-	G	N
Jap Quail	-	-	-	G	-	-	-	-	-	A	A	M	-	-	H	-	L	-	K	-	Q	-	Y	S	-	G	N
Turkey	-	-	-	G	-	-	-	-	-	A	A	M	-	-	L	-	L	-	N	-	-	-	Y	S	-	G	N
RN Pheasant	G	-	-	G	-	-	-	-	-	A	A	M	-	-	M	-	L	-	N	-	-	-	Y	S	-	G	N
Chachalaca	-	I	-	K	-	-	-	-	-	A	A	M	-	-	Y	-	L	-	N	-	-	-	Y	S	-	G	N
Duck A	-	-	-	-	-	-	-	-	-	A	A	M	-	-	L	-	L	-	N	-	-	-	Y	S	-	G	N
Guinea Fowl	-	-	F	G	-	-	-	-	-	A	A	M	-	-	H	-	L	-	N	-	-	-	Y	S	-	G	N
Genetic Marker Panel 2 (28 Markers)																											
Sample ID	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70
Mouse intest	N	Y	N	R	D	R	S	T	D	Y	G	I	F	Q	I	N	S	R	Y	W	C	N	D	G	K	T	P
Mouse milk	N	Y	N	R	D	Q	S	T	D	Y	G	I	F	Q	I	N	S	R	Y	W	C	N	D	G	K	T	P
HEL	-	R	-	T	-	G	-	-	-	-	-	-	L	-	-	-	-	-	W	-	-	-	-	-	K	R	-
Cal Quail	-	R	-	T	-	G	-	-	-	-	-	V	L	-	-	-	-	-	W	-	-	-	-	-	R	-	-
BW Quail	-	R	-	T	-	G	-	-	-	-	-	V	L	-	-	-	-	-	W	-	-	-	-	-	-	-	-
Jap Quail	-	R	-	T	-	G	-	-	-	-	-	-	L	-	-	-	-	-	W	-	-	-	-	-	R	-	-
Turkey	-	R	-	T	-	G	-	-	-	-	-	-	L	-	-	-	-	-	W	-	-	-	-	-	R	-	-
RN Pheasant	-	R	-	T	-	G	-	-	-	-	-	-	L	-	-	-	-	-	W	-	-	-	-	-	R	-	-
Chachalaca	-	R	S	N	-	G	-	-	-	-	-	-	E	-	-	-	-	-	W	-	-	-	-	-	R	-	-
Duck A	-	R	-	T	-	G	-	-	-	-	-	V	L	-	-	-	-	-	W	-	-	-	-	-	-	-	-
Guinea Fowl	-	R	-	T	-	G	-	-	-	-	-	-	L	-	-	-	-	-	W	-	-	-	-	-	R	-	-
Genetic Marker Panel 3 (28 Markers)																											
Sample ID	87	88	89	90	91	92	93	94	95	96	97	98	99	100	101	102	103	104	105	106	107	108	109	110	111	112	113
Mouse intest	D	I	T	A	A	I	Q	C	A	K	R	V	V	R	D	P	Q	G	I	R	A	W	V	A	W	R	T
Mouse milk	D	I	T	A	A	I	Q	C	A	K	R	V	V	R	D	P	Q	G	I	R	A	W	V	A	W	R	A
HEL	-	-	-	-	S	V	N	-	-	-	K	I	-	S	-	G	N	-	M	N	-	-	-	-	-	-	N
Cal Quail	-	-	-	-	T	V	N	-	-	-	K	I	-	S	-	G	N	-	M	N	-	-	-	-	-	-	N
BW Quail	-	-	-	-	T	V	N	-	-	-	K	I	-	S	-	G	N	-	M	N	-	-	-	-	-	-	N
Jap Quail	-	-	-	-	S	V	N	-	-	-	K	I	-	S	-	V	H	-	M	N	-	-	-	-	-	-	N
Turkey	-	-	-	-	S	V	N	-	-	-	K	I	A	S	G	G	N	-	M	N	-	-	-	-	-	-	N
RN Pheasant	-	-	-	-	S	V	N	-	-	-	K	I	-	S	-	G	N	-	M	N	-	-	-	-	-	K	H
Chachalaca	-	A	P	-	S	V	R	-	-	-	-	I	-	S	-	G	N	-	M	N	-	-	-	-	-	K	H
Duck A	-	-	E	-	-	V	R	-	-	-	-	I	-	S	-	G	N	-	M	N	-	-	-	-	-	-	N
Guinea Fowl	-	-	-	T	T	A	N	-	-	-	K	I	-	S	-	G	N	-	M	N	-	-	-	-	-	K	H

3D molecular analysis of the paratope-epitope interface

A “Structural Epitope” consists
of amino acid residues that
contact all six CDRs of antibody

It involves 15-25 residues

Which Residues are Critical for Antigen-Antibody Contact?

- Antibody reactivity with natural variants
- Effect of residue substitutions on antigen-antibody binding
 - Site mutagenesis (alanine scanning)
 - Identification of “energetic” contact residues (“hot spots”)
 - Which energetic residues are non-self?

This analysis was done on data published in the literature


Two Types of Residues in a Structural Epitope

Most residues of a structural epitope can be replaced by other residues without impairing its binding with antibody


Some residues cannot be replaced without causing a major loss in binding activity, they constitute the “hot spot” or “functional” part of the epitope.

Energetic Residues of Functional HEL-D1.3 Epitope, the Corresponding Functional Paratope and in the Epitope-Paratope Interface


HEL-D1.3 Ep



HyHEL-D1.3 CDR




Interface




119D, 124I and 125R are within 3.0 Angstroms from 121Q^{NS}

Functional HEL-5 Epitope


HEL-5 Ep



HyHEL-5 CDR




Interface




45R^{NS} and 68R^{NS} are
3.5 Angstroms apart


HEL-63 Ep



HyHEL-63 CDR




Interface





20Y is 3.5 Angstroms
away from 97K^{NS}

INF γ R-A6 Ep



INF γ R-A6 CDR



All functional epitope residues are within
3.0-3.5 Angstroms from each other


These functional epitopes (“hot spots”) are represented by patches of few residues that are 3.0-3.5 Angstroms apart

Residues are in linear or discontinuous sequences


Other epitopes are represented
by two distinct patches

Examples of HEL Epitopes Defined by a Pair of Patches


HEL-F9.13.7 Ep



HEL-F9.13.7 CDRs



Interface




Patch 1: 73R^{ns}



Patch 2: 20Y and 96K^{ns}

Distance: 13 Angstroms

HyHEL-10 Ep



HyHEL-10 CDRs




Patch 1: 93N^{NS} and 97K^{NS}

Patch 2: 20Y and 21R (self)


Distance: 7.5 Angstroms

Two-Patch Shape of an Epitope on Horse Cytochrome C Recognized by a Mouse Monoclonal Antibody


hCytC-E8 Ep



hCYTC-E8 CDR



Interface



Patch 1: 60K^{NS} and 62E^{NS}
Patch 2: 103N and 104E (self)
Patches are 11.5 Angstroms apart

Only Six Non-Self Residues on Horse Cytochrome C

Horse:

2 GDVEKGKKIFVQKCAQCHTVEKGGKHKTGPNLHGLFGRKTGQAPGF**TY**TDANKNKGITW**K** 61

Mouse:

2 GDVEKGKKIFVQKCAQCHTVEKGGKHKTGPNLHGLFGRKTGQAAGFSYTDANKNKGITWG 61

Horse

62 E**E**TLMEYLENPKKYIPGTKMIFAGIKKK**T**ER**E**DLIAYLK**K**ATNE 105


Mouse

62 EDTLMEYLENPKKYIPGTKMIFAGIKKKGERADLIAYLK**K**ATNE 105

How Do Antibodies Bind to Epitopes?


- So-called “**structural**” epitopes on protein antigens have 15-25 surface residues that can make contact with all six CDRs of antibody
- These residues cover an area of 700-900 square Ångstroms
- Within each lies a centrally located “**functional**” epitope with a few residues and which binds to CDR-H3 plays a dominant role in specific recognition

CDRs of Antibody



- Complementarity Determining Regions H1, H2, H3, on heavy chain and L1, L2 and L3 on light chain represent the “binding face” or “paratope” of an antibody
- The antibody specificity is often determined by a centrally located loop (CDR-H3) that binds to a short amino acid sequence of the protein antigen
- Other CDRs serve as contact sites to stabilize binding to antigen (they play a role in the affinity of antibody)


Antigen-Antibody Complex




How Do Antibodies Bind to Epitopes on Proteins?

- So-called “**structural**” epitopes on protein antigens have 15-25 surface residues that can make contact with all six CDRs of antibody
- These residues cover an area of 700-900 square Ångstroms
- Within each lies a centrally located “**functional**” epitope with a few residues and which binds to CDR-H3 plays a dominant role in specific recognition

Antigen-Antibody Complex




CDRs of Antibody



HEL antibodies Hen Egg Lysozyme Epitopes

Antibody
CDR loops
L1: light blue
L2: cyan
L3: dark blue
H1: orange
H2: magenta
H3: red



Antigen
residues
contacted by

Light chains:
orange



Heavy chains:
blue

Both chains:
green

Structural Aspects of Complexes with HLA Molecules

- Crystal models of MHC-TCR complexes
- Crystal structure of antibody-HLA + peptide complex

T-Cell Receptor-HLA Complex



TCR


HLA

Crystal Structure of HLA-A1-MAGE-A1 Complex with Antibody Fab-Hyb3


- Structural epitope: 4 peptide residues + 9 HLA-A1 residues in α helices
- Paratope: 11 VH and 4 VL residues
- Predominant role of CDR-H2 and H3 in peptide contact

Molecular modeling of the crystalline structure of human monoclonal antibody Hyb3 complexed with HLA-A1-MAGE-A1 as reported by Ziegler's group*


Antigen-Antibody Complex



Structural paratope of Hyb3



Structural epitope of HLA-A1-MAGE-A1



Contact residues are within a molecular surface of about 900 Å²

*Hulsmeyer et al. J Biol. Chem. 280:2972-2980, 2005

Crystal Structure of HLA-A1-MAGE-A1 Complex with Antibody Fab-Hyb3


- Structural epitope: 4 peptide residues + 9 HLA-A1 residues in α helices
- Paratope: 11 VH and 4 VL residues
- Predominant role of CDR-H2 and H3 in peptide contact
- Functional epitope:
 - Contact sites: 65R and 72Q of HLA-A1
 - Antibody does not react with MAGE-A1 presented by HLA-A3 (only difference is contact residue 76A)
 - Antibody does not react with MAGE-A3 which has three residue differences
 - No mutagenesis information

Hulsmeyer et al J. Biol. Chem. 280: 2972-2790, 2005
(Andrew Ziegler's lab at Humboldt University, Berlin)


Prominent residues of the structural MAGE-A1-HLA-A1-Hyb3 epitope and their interface with the Hyb3 paratope

(adapted from Hulsmeyer et al, 2005)

epitope



paratope






Patches are 6-8 Angstroms away

How do we describe HLA epitopes?




- Public and private HLA determinants
(Cross-Reacting Groups, CREG)

Polymorphic Residues on Class I Antigen




HLA-A2



HLA-B27



HLA-Cw3



How do we describe HLA epitopes?

- Public and private HLA determinants (Cross-Reacting Groups, CREG)
- Single amino acid residues
- HLAMatchmaker: Triplets (2002)
- HLAMatchmaker: Eplets (2006)

Concepts about Antibody-Antigen Interactions

- Contact areas of CDR loops of antibody and epitopes on protein antigens (paratope-epitope interface)
- Determination of “structural” epitopes
- Identification of “functional” epitopes
- Determine the sizes and shapes of functional epitopes

Apply concepts to eplet version of HLAMatchmaker

Functional epitopes are
represented by patches of few
residues that are 3.0-3.5
Angstroms apart

In HLAMatchmaker:
Eplets are patches consisting of a
polymorphic residue and residues
within a 3.0-3.5 Angstrom radius

Recent Update:
Designing Epitope Repertoires for
HLAMatchmaker and the HLA
Epitope Registry

Reference:

Update of the HLA Class I Eplet Database in the
Website Based Registry of Antibody-Defined HLA
Epitopes, Tissue Antigens, 83:382-390, 2014

Eplets

- Essential components of HLA epitopes recognized by antibody
- Amino acid configurations within a 3 Angstrom radius of surface-exposed polymorphic residues
- Parts of “structural” HLA epitopes that contact the CDRs of antibody

Design

1. With Cn3D modelling of HLA molecules identify all combinations of polymorphic residues within a 3 Ångstrom radius

There are 116 configurations for exposed residues and 21 configurations solely defined by hidden polymorphic residues

exposed	seq	seq	seq	seq	exposed	seq	seq	seq	seq	exposed	seq	seq	seq	seq	hidden	seq	seq
1	1				73	73	74	76		149	149	151			4	4	
14	14	17			73	73	76			150	150	151	152		6	6	
16	16				73	73	76	77		150	150	151			9	9	
19	19				73	73				150	150				11	11	12
41	41				76	76	77			151	151				21	21	
44	44	45	46		76	76	77	79		151	151	152			24	24	
45	45	46			76	76	79			156	156				30	30	
56	56				76	76	77	80		156	156	158			32	32	
62	62	63			76	76	80			158	158				35	35	
62	62	63	65		76	76	79	80		161	161				49	49	
62	62	63	66		76	76				162	162	163	167		67	67	
62	62	65			77	77	79	80		163	163				70	70	
62	62	65	66		80	80	81	82		163	163	166			74	74	
62	62	66			80	80				163	163	166	167		94	94	95
62	62				80	80	81			166	166	167			97	97	
63	63	65			80	80	82	83		166	166				99	99	
63	63	66			82	81	82	83		170	170	171			113	113	114
63	63	65	66		82	82	83			173	173				116	116	
65	65	66			82	82				177	177				152	152	
65	65	66	67		90	90	91			177	177	178			156	156	
65	65	66	69		102	102	103			180	180				199	199	
65	65	69			105	105				182	182	183	184				
65	65				107	107				186	186						
66	66	67			109	109				193	193	194					
66	66	67	69		127	127				193	193						
66	66	69			131	131				194	194						
66	66	69	70		138	138				207	207						
66	66				138	138	142			211	211						
69	69	70	71		142	142	143	144		219	219						
69	69	70			142	142	145			245	245	246					
69	69	70	73		144	144	145			247	247						
69	69	71			144	144	149			248	248						
69	69	73			144	144				249	249						
69	69	73	77		145	145	149			253	253						
69	69				145	145				261	261	270					
70	70	73			147	147				211	211	212					
70	70				149	149	150	151		273	273						
71	71	73	77		149	149	150			275	275	276					
73	73	74			149	149											

**Configurations of
polymorphic
sequence positions on
class I HLA molecules**

Design

1. With Cn3D modelling of HLA molecules identify all combinations of polymorphic residues within a 3 Ångstrom radius

There are 116 configurations for exposed residues and 21 configurations solely defined by hidden polymorphic residues

2. With HLA Patch Generator (developed by **Dr. Grzegorz Dudek** (now at Czestchowa University of Technology, Poland).

determine polymorphic residue compositions of 3 Ångstrom patches from amino acid sequences of 107 HLA alleles in commercially available Luminex panels

Results: 616 patches were identified, 82 of them are monomorphic for one or two ABC loci and they were ruled out as potential epitopes. Many of the remaining 534 patches are shared by identical groups of Luminex alleles.

Multiple patches shared between the same group of Luminex alleles

Patch	Luminex Alleles						
73TAN	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
73TA	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76ANG	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AGT	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AG	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AT	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76A	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AN	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76ANT	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01


Patches shared between the same group of Luminex alleles

Patch	Luminex Alleles						
73TAN	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
73TA	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76ANG	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AGT	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AG	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AT	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76A	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76AN	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01
76ANT	A*01:01	A*26:01	A*29:01	A*29:02	A*36:01	A*43:01	A*80:01

Eplet: 76ANT Residue description: 73T76A77N79G80T


Patches shared between A*01:01 and A*36:01

Patch	Luminex Alleles	
44KM	A*01:01	A*36:01
149AVH	A*01:01	A*36:01
149AV	A*01:01	A*36:01
150V	A*01:01	A*36:01
150VH	A*01:01	A*36:01
150VHA	A*01:01	A*36:01
156RV	A*01:01	A*36:01
158V	A*01:01	A*36:01



Patches shared between A*01:01 and A*36:01

Patch	Luminex Alleles	
44KM	A*01:01	A*36:01
149AVH	A*01:01	A*36:01
149AV	A*01:01	A*36:01
150V	A*01:01	A*36:01
150VH	A*01:01	A*36:01
150VHA	A*01:01	A*36:01
156RV	A*01:01	A*36:01
158V	A*01:01	A*36:01



Eplet: 44KM₃ Description: 44K45M(149A150V151H152A)(156R158V)

Overlapping eplets that are shared between similar but not identical groups of Luminex alleles

Example: six 62R-related eplets on Luminex alleles

[illegible]

Design

1. With Cn3D modelling of HLA molecules identify all combinations of polymorphic residues within a 3 Ångstrom radius

There are 116 configurations for exposed residues and 21 configurations solely defined by hidden polymorphic residues

2. With HLA Patch Generator (developed by Dr. Grzegorz Dudek, Czestchowa University of Technology, Poland) determine polymorphic residue compositions of 3 Ångstrom patches from amino acid sequences of 107 HLA alleles in commercially available Luminex panels

Results: 616 patches were identified, 82 of them are monomorphic for one or two class I loci and were ruled out as potential epitopes. Many of the remaining 534 patches are shared by identical groups of Luminex alleles

3. Annotations of 219 eplets defined by polymorphic residues on the molecular surface and 51 eplets defined solely by polymorphic residues in antibody-inaccessible positions

219 Eplets in antibody-accessible positions on the molecular surface of HLA-A, B, C antigens


Eplet	#A	#B	#C	Total	Eplet	#A	#B	#C	Total	Eplet	#A	#B	#C	Total	Eplet	#A	#B	#C	Total	Eplet	#A	#B	#C	Total	Eplet	#A	#B	#C	Total
1C	0	0	17	17	71ATN	0	4	4	8	145KHA	7	0	0	7	63NN	10	0	0	10	76VDT	20	0	0	20	163T	25	16	17	58
16S	0	0	3	3	71ATS	0	11	8	19	145RT	7	0	0	7	65GK	4	0	0	4	76VRN	0	2	13	15	163TEW	22	16	17	55
17RS	2	0	0	2	71HS	22	0	0	22	147L	0	3	4	7	65QIA	0	13	0	13	76VS	0	1	13	14	163TG	3	0	0	3
17WR	0	0	1	1	71KA	0	4	0	4	149AH	14	0	0	14	65QKR	0	1	21	22	77NGT	8	0	0	8	166DG	5	1	0	6
19K	1	0	0	1	71QS	13	0	0	13	149TAH	8	0	0	8	65QNR	0	0	2	2	77SRN	0	38	13	51	166ES	0	5	0	5
41T	0	13	0	13	71SA	0	4	0	4	150AAH	12	0	0	12	65RA	30	4	0	34	79GT	28	0	0	28	170RH	1	9	0	10
43RRM	2	0	0	2	71STN	12	0	0	12	150AH	20	0	0	20	65RK	6	0	0	6	80I	6	12	0	18	170RY	35	50	16	99
44KM ₃	2	0	0	2	71TTD	0	2	0	2	151AHA	2	0	0	2	65RNA	24	4	0	28	80K	0	0	10	10	173K	0	0	3	3
44RM	33	12	0	45	71TTN	0	12	0	12	151AHE	9	0	0	9	66I	0	54	0	54	80N	0	39	13	52	177DK	0	5	0	5
44RMA	0	11	0	11	71TTS	0	27	0	27	151AHV	9	0	0	9	66IC	0	13	0	13	80T	28	8	0	36	177DT	0	3	0	3
44RME	33	0	0	33	72QTD	25	22	12	59	151ARV	9	35	0	44	66ICT	0	9	0	9	80TA	0	4	0	4	177KT	0	0	4	4
44RT	0	10	0	10	73ID	3	0	0	3	151H	22	0	0	22	66IF	0	8	0	8	80TL	28	4	0	32	180E	0	8	0	8
45EE	0	26	0	26	73AN	0	0	6	6	152HA	4	0	0	4	66IS	0	22	0	22	80TLR	0	8	0	8	182TDP	15	0	0	15
45KE	0	11	0	11	73AS	0	0	5	5	152RA	0	0	4	4	66IY	0	11	0	11	81ALR	6	16	0	22	184A	19	0	0	19
56E ₄	1	0	0	1	73TD	17	4	0	21	152RE	0	24	18	42	66IYT	0	2	0	2	82LR	6	20	0	26	184H	0	0	19	19
56R	3	0	0	3	73TDA	6	0	0	6	152RR	2	0	0	2	66K	10	1	21	32	90AR	0	0	1	1	184R ₃	0	0	1	1
62EE	5	0	0	5	73TDE	7	20	0	27	152RT	0	0	1	1	66KA	10	0	0	10	90D	11	1	8	20	186R	1	0	0	1
62GE	5	3	0	8	73TDV	12	2	12	26	152RW	1	0	0	1	66KAH	9	0	0	9	102DV	33	42	3	78	193AV	23	0	0	23
62GK ₂	5	0	0	5	73TN	12	16	4	32	156DA	0	8	1	9	66N	24	4	2	30	102HV	1	0	0	1	193LV	0	0	1	1
62GRN	0	3	0	3	73TS	2	37	9	48	156LA	16	40	8	64	66NAH	12	0	0	12	103L	0	16	20	36	193PI	12	51	0	63
62LQ	3	0	0	3	73TV	17	2	12	31	156QA	4	0	1	5	66NAQ	12	0	0	12	103M ₂	0	1	0	1	193PL ₃	0	0	3	3
62QE	10	0	0	10	73TVD	17	0	0	17	156RA	0	3	8	11	66NM	2	4	0	6	105S	21	0	0	21	193PV	0	8	19	27
62REN	0	1	2	3	73TVN	0	0	4	4	156WA	12	4	5	21	66NV	23	0	0	23	107W	6	0	0	6	194V	22	8	20	50
62RER	0	1	0	1	73TVS	0	1	8	9	158T	0	4	0	4	66RKQ	1	0	0	1	109F	32	0	0	32	207S	23	0	0	23
62RK	1	1	21	23	73TY	0	37	0	37	161D	1	0	0	1	69AA	0	17	0	17	109FE	31	0	0	31	211T	0	0	1	1
62RN	11	32	0	43	76ANT	7	0	0	7	162DLS	0	1	0	1	69AQT	13	9	0	22	127K	12	0	0	12	219W	0	0	9	9
62RNQ	0	32	0	32	76ED	0	4	0	4	162GLS	0	4	0	4	69AT	31	17	0	48	131S	0	47	0	47	245AS	14	0	0	14
62RNR	11	0	0	11	76EG	1	0	0	1	163E	2	14	3	19	69ATD	17	2	0	19	138K	0	0	2	2	245TA	0	2	0	2
62RR	11	1	0	12	76EN	5	16	0	21	163EW	1	14	3	18	69ATN	12	4	0	16	138MI	26	0	0	26	245VA	2	0	0	2
62RRN	10	1	0	11	76ENI	4	12	0	16	163L	0	29	3	32	69ATS	2	10	0	12	142ITQ	18	56	22	96	248M	0	0	1	1
62RTN	10	1	2	13	76ENR	4	16	0	20	163LE	0	28	3	31	69RA	0	0	11	11	143S	0	3	1	4	253Q	22	1	4	27
63EI	0	22	0	22	76ENT	1	4	0	5	163LG	0	1	0	1	69RT	0	1	12	13	144K	16	0	0	16	267PE	1	0	19	20
63EK	9	1	21	31	76ES	2	37	0	39	163LW	0	23	3	26	69TNT	0	41	0	41	144KA	15	0	0	15	267QE	0	1	4	5
63EN	11	4	2	17	76ESI	2	0	0	2	163R	7	0	0	7	70HT	18	0	0	18	144QL	0	2	0	2	270C	0	1	1	2
63ER	16	4	0	20	76ESN	0	37	0	37	163RG	1	0	0	1	70IAQ	0	9	0	9	144TKH	8	0	0	8	275EL	7	0	0	7
63ERN	12	4	0	16	76ET	1	8	0	9	163RW	6	0	0	6	70QT	13	10	12	35	144KR	8	0	0	8	275G	0	0	3	3
															71ATD	0	2	0	2	145HT	1	0	0	1	275K	0	1	5	6

51 Eplets based solely on residue polymorphisms in antibody-inaccessible sequence positions

<u>Eplet</u>	<u>#A</u>	<u>#B</u>	<u>#C</u>	<u>Total</u>	<u>Eplet</u>	<u>#A</u>	<u>#B</u>	<u>#C</u>	<u>Total</u>	<u>Eplet</u>	<u>#A</u>	<u>#B</u>	<u>#C</u>	<u>Total</u>
9D	0	1	6	7	94IL	0	0	1	1	113HD	0	16	1	17
9F	10	0	1	11	94TF	0	0	1	1	113HN	0	26	0	26
9H	0	14	0	14	94TI	25	1	1	27	113YD	0	4	14	18
9S	6	0	2	8	94TL	6	31	17	54	113YH	11	4	0	15
9T	5	0	0	5	94TV	4	0	0	4	113YN	0	9	8	17
9Y	14	44	14	72	94TW	0	16	0	16	113YQ	12	0	0	12
11AM	0	37	0	37	97I	9	0	0	9	113YR	10	0	0	10
11AV	0	3	20	23	97M	13	0	0	13	116D	22	6	0	28
11SM	1	0	0	1	97N	0	3	0	3	116F	0	10	9	19
11SV	34	19	3	56	97R	13	32	18	63	116L	0	11	1	12
21H	0	0	7	7	97S	0	9	0	9	116S	0	15	10	25
24S	0	21	7	28	97T	0	11	0	11	116Y	11	17	3	31
24T	0	17	0	17	97V	0	2	0	2	152A	4	0	4	8
30G	0	1	0	1	97W	0	2	5	7	152E	9	24	18	51
32L	0	14	0	14	99F	4	3	5	12	152V	19	35	0	54
35Q	1	0	3	4	99S	0	1	1	2	156L	17	44	8	69
94II	0	11	3	14	99Y	31	55	13	99	156R	2	3	8	13
										199V	0	2	0	2

Eplets on Non-Luminex Alleles

Comparisons of relevant amino acid configurations between 76ANT-carrying Luminex alleles and non-Luminex alleles




Selected Alleles		73	76	77	79	80	Selected Alleles		73	76	77	79	80	Selected Alleles		73	76	77	79	80
A*01:01	LUM	T	A	N	G	T	A*26:06	T	V	D	G	T		A*26:37		T	A	N	G	T
A*01:02		T	A	N	G	T	A*26:07	T	A	N	G	T		A*29:01	LUM	T	A	N	G	T
A*01:03		T	A	N	G	T	A*26:08	T	A	N	G	T		A*29:02	LUM	T	A	N	G	T
A*01:06		T	A	N	G	T	A*26:09	T	A	N	G	T		A*29:03		T	A	N	G	T
A*01:07		T	E	N	G	T	A*26:10	T	A	N	G	T		A*29:04		T	A	N	G	T
A*01:08		T	A	N	G	T	A*26:12	T	A	N	G	T		A*29:05		T	A	N	G	T
A*01:09		T	A	N	G	T	A*26:13	T	A	N	G	T		A*29:06		T	A	N	G	T
A*01:10		T	A	N	G	T	A*26:14	T	A	N	G	T		A*29:07		T	A	N	G	T
A*01:12		T	A	N	G	T	A*26:15	T	A	N	G	T		A*29:09		T	A	N	G	T
A*01:13		T	V	D	G	T	A*26:16	T	A	N	G	T		A*29:10		T	A	N	G	T
A*01:14		T	A	N	G	T	A*26:17	T	A	N	G	T		A*29:11		T	A	N	G	T
A*01:17		T	A	N	G	T	A*26:18	T	A	N	G	T		A*29:12		T	A	N	G	T
A*01:19		T	A	N	G	T	A*26:19	T	A	N	G	T		A*29:13		T	E	S	R	I
A*01:20		T	A	N	G	T	A*26:20	T	A	N	G	T		A*29:14		T	A	N	G	T
A*01:21		T	A	N	G	T	A*26:21	T	V	D	G	T		A*29:15		T	A	N	G	T
A*01:23		T	A	N	G	T	A*26:22	T	A	N	G	T		A*29:16		T	A	N	G	T
A*01:24		T	A	N	G	T	A*26:23	T	A	N	G	T		A*29:17		T	A	N	G	T
A*01:25		T	A	N	G	T	A*26:24	T	A	N	G	T		A*29:18		T	A	N	G	T
A*01:26		T	A	N	G	T	A*26:26	T	A	N	G	T		A*29:19		T	V	D	G	T
A*01:28		T	V	D	G	T	A*26:27	T	A	N	G	T		A*36:01	LUM	T	A	N	G	T
A*01:29		T	A	N	G	T	A*26:28	T	A	N	G	T		A*36:02		T	A	N	G	T
A*01:30		T	A	N	G	T	A*26:29	T	A	N	G	T		A*36:03		T	A	N	G	T
A*01:32		T	A	N	G	T	A*26:30	T	V	D	G	T		A*36:04		T	A	N	G	T
A*01:33		T	A	N	G	T	A*26:31	T	A	N	G	T		A*43:01	LUM	T	A	N	G	T
A*26:01	LUM	T	A	N	G	T	A*26:32	T	A	N	G	T		A*80:01	LUM	T	A	N	G	T
A*26:02		T	A	N	G	T	A*26:33	I	A	N	G	T		A*11:17		T	A	N	G	T
A*26:03		T	V	D	G	T	A*26:34	T	A	N	G	T		A*11:40		T	A	N	G	T
A*26:04		T	A	N	G	T	A*26:35	T	A	N	G	T		A*24:04		T	A	N	G	T
A*26:05		T	E	N	G	T	A*26:36	T	A	N	G	T		A*74:10		T	A	N	G	T

Residue configurations of non-Luminex alleles in relation to an eplet annotated as 44KM₃ and shared between A*01:01 and A*36:01

Alleles			44	45	149	150	151	152	156	158
A*01:01 LUM	44KM ₃		K	M	A	V	H	A	R	V
A*36:01 LUM	44KM ₃		K	M	A	V	H	A	R	V
A*01:02	44KM ₃		K	M	A	V	H	A	R	V
A*01:03	44KM ₃		K	M	A	V	H	A	R	V
A*01:06			K	M	A	V	H	A	L	A
A*01:07	44KM ₃		K	M	A	V	H	A	R	V
A*01:08	44KM ₃		K	M	A	V	H	A	R	V
A*01:09	44KM ₃		K	M	A	V	H	A	R	V
A*01:10			K	M	A	A	R	R	R	V
A*01:12			K	M	A	A	H	V	Q	A
A*01:13	44KM ₃		K	M	A	V	H	A	R	V
A*01:14	44KM ₃		K	M	A	V	H	A	R	V
A*01:17	44KM ₃		K	M	A	V	H	A	R	V
A*01:19			K	M	A	A	H	V	Q	A
A*01:20	44KM ₃		K	M	A	V	H	A	R	V
A*01:21			K	M	A	A	H	V	R	V
A*01:23	44KM ₃		K	M	A	V	H	A	R	V
A*01:24	44KM ₃		K	M	A	V	H	A	R	V
A*01:25			K	M	A	V	H	A	Q	V
A*01:26			K	M	A	A	H	A	R	V
A*01:28	44KM ₃		K	M	A	V	H	A	R	V
A*01:29	44KM ₃		K	M	A	V	H	A	R	V
A*01:30	44KM ₃		K	M	A	V	H	A	R	V
A*01:32	44KM ₃		K	M	A	V	H	A	R	V
A*01:33	44KM ₃		K	M	A	V	H	A	R	V
A*01:35	44KM ₃		K	M	A	V	H	A	R	V
A*36:02			K	M	A	V	H	A	R	A
A*36:03	44KM ₃		K	M	A	V	H	A	R	V
A*36:04	44KM ₃		K	M	A	V	H	A	R	V
A*03:18			R	M	A	V	H	A	R	V
A*68:41			R	M	A	V	H	V	W	A

Residue configurations of non-Luminex alleles in relation to the 76ED and 80TLL eplets shared between Luminex alleles B*27:03, B*27:05, B*37:01 and B*47:01

<u>Alleles</u>	<u>76ED/80TLL Eplet</u>	<u>76</u>	<u>77</u>	<u>80</u>	<u>81</u>	<u>82</u>
B*27:03 LUM	76ED and/or 80TLL	E	D	T	L	L
B*27:05 LUM	76ED and/or 80TLL	E	D	T	L	L
B*37:01 LUM	76ED and/or 80TLL	E	D	T	L	L
B*47:01 LUM	76ED and/or 80TLL	E	D	T	L	L
B*27:07/09/10/13/14/16/17/19/27/28/ 29/32/34/35/37/38/39/41/43/45/46/50	76ED and/or 80TLL	E	D	T	L	L
B*37:02/04/06/07/08/09/12/13, B*47:05	76ED and/or 80TLL	E	D	T	L	L
B*07:27, B*15:43, B*38:17, B*53:03	76ED and/or 80TLL	E	D	T	L	L
B*27:42, B*37:05	76ED only	E	D	N	L	R
B*27:04/06/11/15/20/21/23/24/25/31/36	80TLL only	E	S	T	L	L
B*27:08 LUM	not 76ED/80TLL	E	S	N	L	R
B*27:12/18/26/33/40/44, B*37:11, B*47:02	not 76ED/80TLL	E	S	N	L	R
B*37:34, B*47:03	not 76ED/80TLL	E	S	N	L	L
B*27:01/02/30, B*37:10, B*47:04	not 76ED/80TLL	E	N	T	A	L



Distance between 76D and 82L is 7 Ångstroms


HLAMatchmaker:

Examples of Eplets in Class I Sequence Positions





	44	62	65	69	77	...	163	184	193	248
A*0101	44KM	62QE	65RNM	66RNH	77ANT		163RG	184P	193PI	248VK
A*0201	44RM	62GE	65RKV	70KAH	77VDI		163TW	184A	193AV	248VK
A*0202	44RM	62GE	65RKV	70KAH	77VDI		163TW	184A	193AV	248VK
A*0203	44RM	62GE	65RKV	70KAH	77VDI		163TW	184A	193AV	248VK
A*0205	44RM	62GE	65RKV	70KAH	77VDI		163TW	184A	193AV	248VK
A*0301	44RM	62QE	65RNV	66RNQ	77VDI		163TW	184P	193PI	248VK
A*0302	44RM	62QE	65RNV	66RNQ	77VDI		163TW	184P	193PI	248VK
A*1101	44RM	62QE	65RNV	66RNQ	77VDI		163RW	184P	193PI	248VK
A*1102	44RM	62QE	65RNV	66RNQ	77VDI		163RW	184P	193PI	248VK
A*2301	44RM	62EE	66GKH	70KAH	77ENI		163TG	184P	193PI	248VK
A*2402	44RM	62EE	66GKH	70KAH	77ENI		163TG	184P	193PI	248VK
A*2403	44RM	62EE	66GKH	70KAH	77ENI		163TW	184P	193PI	248VK
A*2501	44RM	62RN	65RNV	66RNH	77ESI		163RW	184A	193AV	248VK
A*2601	44RM	62RN	65RNV	66RNH	77ANT		163RW	184A	193AV	248VK
A*2901	44RM	62LQ	65RNV	66RNQ	77ANT		163TW	184A	193AV	248VK
A*2902	44RM	62LQ	65RNV	66RNQ	77ANT		163TW	184A	193AV	248VK
A*3001	44RM	62QE	65RNV	66RNQ	77VDI		163TW	184P	193PI	248VK
A*3002	44RM	62QE	65RNV	66RNH	77ENT		163TW	184P	193PI	248VK
A*3101	44RM	62QE	65RNV	66RNH	77VDI		163TW	184P	193AV	248VK

HLA Epitopes

- HLAMatchmaker considers **eplets** as equivalents to **functional** epitopes
- Amino acid residues within a 15 Ångstrom radius of eplets can contribute **structural** HLA epitopes



Four Examples on
HLA-A2 molecules



Two strategies to determine the antibody-verified HLA epitope repertoire

1. HLAMatchmaker is a **theoretical algorithm that predicts** HLA epitopes on the HLA molecular surface from stereochemical modeling of epitope-paratope interfaces of antigen-antibody complexes
2. An **empirical method** to analyze antibody reactivity with HLA-typed panels and the identification of amino acid configurations shared between reactive alleles

Two strategies to determine the HLA epitope repertoire

1. HLAMatchmaker is a **theoretical algorithm that predicts** HLA epitopes on the HLA molecular surface from stereochemical modeling of epitope-paratope interfaces of antigen-antibody complexes
2. An **empirical method** to analyze antibody reactivity with HLA-typed panels and the identification of amino acid configurations shared between reactive alleles

Terasaki and his group used single allele Luminex panels tested with mouse monoclonal antibodies and (absorbed/eluted) alloantibodies
They reported 103 HLA-A,B,C epitopes annotated with a numbering system
(El-Awar et al. : *Transplantation* 84:532-540, 2007)

How do “Terasaki’s” epitopes correspond to eplets ?

Tissue Antigens, 74:117-133, 2009; 74:134-146, 2009

Results

50/103 Terasaki's class I
epitopes correspond to eplets


Examples of Terasaki's epitopes that correspond to eplets ^a

<u>TerEp</u>	<u>Antibody Reactive Antigens</u>	<u>Residue Description of TerEp</u>	<u>Eplet</u>
#19	A2,23,24,68,69	127K	127K
#3	A23, 24	65G	65GKA
#28	A23,24,80	62E	62EE
#14	A1,23,24,02,80; B76	166D/ 167G	167DG
#24	A23,24,25,32; B13,2705,37,38, 44,47,49,51,52,53,57,58,59,63,77	82L/ 83R	82LR
#23	A23,24,25,32; B38,49,51, 52,53,57,58,59,63,77	80I	79RI
#217	B13,2705,37,44,47	76E+80T/ 79R+80T/ 80T+82L/ 80T+83R	79RT
#245	B35,4005,46,49,50,51,52,53,56, 57,58,62,63,71,72,75, 77,78; Cw9,w10	163L+167W	163LW
#35	B18,35,37,51,52,53,58,78	45T	44RT
#222	A6602; B7,13,27,47,48,60,61,73,81; Cw2,w17	163E+166E/ 163E+167W	163EW
#39	Cw2,w9,w10,w15	21H	21H
#246	B46,73; Cw1,w7,w8,w9,w10,w12,w14,w16	76V+80N/ 73T+76V+79R	80VRN
#421	B46; Cw1,w8,w9,w10,w14,w16	(73T)+76V+80N+90A	77TVS


^a These epitopes are on A24, B57 and Cw3

Molecular locations of Terasaki's epitopes and corresponding eplets


A*2402



B*5701



Cw*0303



Epitopes Defined by Eplet Pairs

- Duquesnoy, Mulder, et al: *HLAMatchmaker-Based Analysis of Human Monoclonal Antibody Reactivity Demonstrates the Importance of an Additional Contact Site for Specific Recognition of Triplet-Defined Epitopes* *Human Immunology* 66: 749-761 2005

Antibody Producer: A2,A68; B7,B27; Cw2,Cw7
Mismatched triplets on immunizing A3:
62QE,142MI,144TKR,151AHE,163DT

Hu mAb	Specific epitope
OK2H12	62QE + self 56G
OK4F9 and OKF10	142MI + self 79GTLRG
OK5A3	144TKR + self 151H

Reactivity of mAb OK2H12 with 62Qe-carrying alleles


Allele	Triplets shared with HLA-A3	N	Lymphocytotoxicity		Antigen-Binding	
			%Pos Rx	Rx Strength	Flow Beads	Elisa
A*0301	62Qe,142ml,144tKr,151aHe,163dT	103	99%	7.6	++	++
A*0101	62Qe,142ml,144tKr	61	100%	7.6	++	+
A*0103	62Qe,142ml,144tKr	4	75%	6.3	nd	nd
A*1101	62Qe,142ml,144tKr	83	96%	7.7	++	++
A*3201	62Qe,142ml	29	100%	7.8	++	++
A*3601	62Qe,142ml,144tKr	6	83%	6.5	nd	+
A*7401	62Qe,142ml	8	100%	7.5	nd	++
A*3001	62Qe,142ml	12	0%	1.0	Neg	Neg
A*3002	62Qe,142ml	16	0%	1.0	nd	nd
A*3101	62Qe,142ml	38	8%	1.6	Neg	Neg
A*2301	142ml	31	0%	1.0	Neg	Neg
A*2402	142ml,144tKr	107	4%	1.2	nd	nd
A*2403	142ml,144tKr	4	0%	1.0	Neg	nd
A*2407	142ml,144tKr	6	0%	1.0	nd	nd
A*2501	142ml,151aHe	11	0%	1.0	Neg	Neg
A*2601	142ml,151aHe	36	6%	1.3	Neg	Neg
A*2901	142ml	24	4%	1.2	Neg	Neg
A*3301	142ml	15	0%	1.0	Neg	Neg
A*3303	142ml	16	0%	1.0	nd	Neg
A*3401	142ml,151aHe	22	3%	1.2	nd	Neg
A*3402	142ml,151aHe	11	9%	1.7	Neg	nd
A*6601	142ml,151aHe	14	0%	1.0	Neg	Neg
A*8001	142ml,144tKr	6	0%	1.0	Neg	Neg
A2/A28	none	67	1%	1.1	Neg	Neg

56G is the only residue shared between HLA-A3 and 62Qe-carrying alleles that react with OK2H12


	Rx	9	17	35	44	45	<u>56</u>	<u>62</u>	<u>63</u>	65	66	67	70	73	74	76	77	79	80	81	82	83	90	95	97	99	105	107
A*0301	Pos	F	R	R	R	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	V	Q	T	D	V	D	G	T	L	R	G	A	I	I	Y	S	G
A*0101	Pos	F	R	R	K	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	M	H	T	D	A	N	G	T	L	R	G	D	I	I	Y	P	G
A*0103	Pos	F	R	R	K	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	M	H	T	D	A	N	G	T	L	R	G	D	I	M	Y	P	G
A*1101	Pos	Y	R	R	R	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	V	Q	T	D	V	D	G	T	L	R	G	D	I	I	Y	P	G
A*3201	Pos	F	R	R	R	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	V	H	T	D	E	S	R	I	A	L	R	A	I	M	Y	P	G
A*3601	Pos	F	R	R	K	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	M	H	T	D	A	N	G	T	L	R	G	D	I	I	Y	P	G
A*7401	Pos	F	R	R	R	M	<u>G</u>	<u>Q</u>	<u>E</u>	R	N	V	H	T	D	V	D	G	T	L	R	G	A	I	M	Y	P	G
A*3001	Neg	S	S	R	R	M	<u>R</u>	<u>Q</u>	<u>E</u>	R	N	V	Q	T	D	V	D	G	T	L	R	G	A	I	I	Y	S	G
A*3002	Neg	S	S	R	R	M	<u>R</u>	<u>Q</u>	<u>E</u>	R	N	V	H	T	D	E	N	G	T	L	R	G	A	I	I	Y	S	G
A*3101	Neg	T	R	R	R	M	<u>R</u>	<u>Q</u>	<u>E</u>	R	N	V	H	I	D	V	D	G	T	L	R	G	A	I	M	Y	S	G

		109	114	116	127	138	142	143	144	145	149	150	151	152	156	158	161	163	166	167	171	184	186	193	194	207	246	253
A*0301	Pos	F	R	D	N	M	I	T	K	R	A	A	H	E	L	A	D	T	E	W	Y	P	K	P	I	G	A	E
A*0101	Pos	F	R	D	N	M	I	T	K	R	A	V	H	A	R	V	E	R	D	G	Y	P	K	P	I	G	A	E
A*0103	Pos	F	R	D	N	M	I	T	K	R	A	V	H	A	R	V	E	R	D	G	Y	P	K	P	I	G	A	E
A*1101	Pos	F	R	D	N	M	I	T	K	R	A	A	H	A	Q	A	E	R	E	W	Y	P	K	P	I	G	A	E
A*3201	Pos	L	Q	D	N	M	I	T	Q	R	A	A	R	V	L	A	E	T	E	W	Y	A	K	A	V	S	S	Q
A*3601	Pos	F	R	D	N	M	I	T	K	R	A	V	H	A	R	V	E	T	E	W	Y	P	K	P	I	G	A	E
A*7401	Pos	L	Q	D	N	M	I	T	Q	R	A	A	R	V	L	A	E	T	E	W	Y	A	K	A	V	S	S	Q
A*3001	Neg	F	E	H	N	M	I	T	Q	R	A	A	R	W	L	A	E	T	E	W	Y	P	K	P	I	G	A	E
A*3002	Neg	F	E	H	N	M	I	T	Q	R	A	A	R	R	L	A	E	T	E	W	Y	P	K	P	I	G	A	E
A*3101	Neg	F	Q	D	N	M	I	T	Q	R	A	A	R	V	L	A	E	T	E	W	Y	P	K	A	V	S	S	Q


Shared Polymorphic Amino Acids in the 62Qe-Defined Specificity Site



Shared Polymorphic Amino Acids in the 62Qe-Defined Specificity Site



Shared Polymorphic Amino Acids in the 62Qe-Defined Specificity Site




Many Residue Substitutions
on HLA do not affect
Reactivity with Alloantibody

Permissive Residue Substitutions that Do Not Inhibit Reactivity of the 62Qe+56G-Specific mAb

	Rx	9	17	35	44	45	56	62	63	65	66	67	70	73	74	76	77	79	80	81	82	83
A*0301	Pos	F	R	R	R	M	<u>G</u>	Q	E	R	N	V	Q	T	D	V	D	G	T	L	R	G
A*0101	Pos	F	R	R	K	M	<u>G</u>	Q	E	R	N	M	H	T	D	A	N	G	T	L	R	G
A*0103	Pos	F	R	R	K	M	<u>G</u>	Q	E	R	N	M	H	T	D	A	N	G	T	L	R	G
A*1101	Pos	Y	R	R	R	M	<u>G</u>	Q	E	R	N	V	Q	T	D	V	D	G	T	L	R	G
A*3201	Pos	F	R	R	R	M	<u>G</u>	Q	E	R	N	V	H	T	D	E	S	R	I	A	L	R
A*3601	Pos	F	R	R	K	M	<u>G</u>	Q	E	R	N	M	H	T	D	A	N	G	T	L	R	G
A*7401	Pos	F	R	R	R	M	<u>G</u>	Q	E	R	N	V	H	T	D	V	D	G	T	L	R	G
A*3001	Neg	S	S	R	R	M	<u>R</u>	Q	E	R	N	V	Q	T	D	V	D	G	T	L	R	G
A*3002	Neg	S	S	R	R	M	<u>R</u>	Q	E	R	N	V	H	T	D	E	N	G	T	L	R	G
A*3101	Neg	T	R	R	R	M	<u>R</u>	Q	E	R	N	V	H	I	D	V	D	G	T	L	R	G
		109	114	116	127	138	142	143	144	145	149	150	151	152	156	158	161	163	166	167	171	184
A*0301	Pos	F	R	D	N	M	I	T	K	R	A	A	H	E	L	A	D	T	E	W	Y	P
A*0101	Pos	F	R	D	N	M	I	T	K	R	A	V	H	A	R	V	E	R	D	G	Y	P
A*0103	Pos	F	R	D	N	M	I	T	K	R	A	V	H	A	R	V	E	R	D	G	Y	P
A*1101	Pos	F	R	D	N	M	I	T	K	R	A	A	H	A	Q	A	E	R	E	W	Y	P
A*3201	Pos	L	Q	D	N	M	I	T	Q	R	A	A	R	V	L	A	E	T	E	W	Y	A
A*3601	Pos	F	R	D	N	M	I	T	K	R	A	V	H	A	R	V	E	T	E	W	Y	P
A*7401	Pos	L	Q	D	N	M	I	T	Q	R	A	A	R	V	L	A	E	T	E	W	Y	A
A*3001	Neg	F	E	H	N	M	I	T	Q	R	A	A	R	W	L	A	E	T	E	W	Y	P
A*3002	Neg	F	E	H	N	M	I	T	Q	R	A	A	R	R	L	A	E	T	E	W	Y	P
A*3101	Neg	F	Q	D	N	M	I	T	Q	R	A	A	R	V	L	A	E	T	E	W	Y	P

Locations of Permissive Polymorphic Residues for 62Qe+s56G-Specific mAb Reactivity



Human Monoclonal Antibody Reactivity with HLA Class I Epitopes Defined by Pairs of Mismatched Eplets and Self Eplets


Rene Duquesnoy¹ Marilyn Marrari¹, Justin Mostecky²,
Arend Mulder³, Frans Claas³ and Ivan Balazs²

¹ University of Pittsburgh Medical Center


²Gen-Probe Transplant Diagnostics, Stamford, CT

³ Leiden University Medical Center, The Netherlands

Monoclonal: HDG11G Antibody Producer: A2, 24; B7,40;Cw7,-
 Immunizer: B35 Mismatched eplets:
 44RT,63NQIN,65QIF,131S,**163LW**,193PV




Monoclonal: HDG11G Antibody Producer: A2, 24; B7,40;Cw7,-
 Immunizer: B35 Mismatched eplets:
 44RT,63NQIN,65QIF,131S,**163LW**,193PV




62QRI is a self eplet on B7 and B40 of the antibody producer

Monoclonal: WK1D12 Antibody producer A1,-; B8,-;Cw7,-
 Immunizer: B27 Mismatched eplets:
 32L,62REI,65QIA,71KA,76EDT,79ERT,82LR,131S,163EW



73TE is a self eplet on B8 of the antibody producer

This antibody reacts with all Bw4 antigens except A25 and B13
Monoclonal VDK8F7 Antibody producer: A3, 31; B35,-;Cw4,-
 Immunizer: B37 Mismatched eplets:
 62REI,66QIN,65QIS,73TTD,76EDT,79ERT,82LR,156DA



145RAA is a self eplet on the alleles of the antibody producer

Marrari M, Mostecky J, Mulder A, Balazs I, Claas F, Duquesnoy R. Human Monoclonal Antibody Reactivity with HLA Class I Epitopes Defined by Pairs of Mismatched Eplets and Self Eplets. Transplantation 2010; 90:1468-72

Antibody	HLA type	Imm	Eplet	Distance	Reactive Antigens
Antibody producer	Specificity				
HDG4B1	A2, 24; B7, 60; Cw7, w10	A32	65RNA+s79ERI	11 Å	A25/A32/B17/B63
MUL9E11	A2, 25; B18, 51; Cw12, 15	B55	69AA-s76E	8 Å	B7/B*27/B42/B54/B55/B56/B57/B58/B63/B67/B81/B82
MUL9F4	A2, 25; B18, 51; Cw12, 15	B55	65QIA-s76ES	7 Å	B7/B*2708/B42/B54/B55/B56/B67/B81/B82
MUL4C8	A2, 25; B18, 51; Cw12, 15	A11	144KR+s151H	7 Å	A1/A3/A11/A24/A36
BRO11F6	A26, 68; B38, 44; Cw7	??	144KR+s151H	7 Å	A1/A3/A11/A24/A36
OK2H12	A2, 28; B7, 27; Cw2, w7	A3	62QE-s56G	11 Å	A1/A3/A11/A32/A36/A74
OK4F10	A2, 28; B7, 27; Cw2, w7	A3	142MI+s79GT	13 Å	A1/A3/A11/A26/A29/A30/A31/A33/A34/A36/A43/A66/A74/A80
OK 4F9	A2, 28; B7, 27; Cw2, w7	A3	142MI+s79GT	13 Å	A1/A3/A11/A26/A29/A30/A31/A33/A34/A36/A43/A66/A74/A80
VDK8C7	A3, 31; B35,-; Cw4,-	B37	82LR-s145RAA	14 Å	All Bw4 except A25 and B13
WK1D12	A1,-; B8,-; Cw7,-	?	163EW+s73TE	16 Å	B7/B13/B27/B47/B48/B40/B81
HA5C2	A1, 32; B8, 61; Cw2, w7	A2?	142TKH-149AAH	3 Å	A2 (not A*0203)/A28
HDG11G12	A2, 24; B7, 60; Cw7, w10	B35	163LW+s62RQI	8 Å	B15/B35/B49/B50/B51/B52/B53/B56/B78
			Average	10 Å	

Summary of Terasaki's epitopes and their corresponding eplets

Locus	TerEps	Eplets	Eplet pairs	Eplets with substitutions	No eplets
ABC ^a	103	50	31	12	10
DR ^b	60	45	3	12	2
DQ ^b	18	13	1	3	1

^a Duquesnoy RJ and Marrari M, *Tissue Antigens*, 74:117-133, **2009**.

^b Marrari M and Duquesnoy RJ, *Tissue Antigens*, 74:134-146, **2009**.

16th International Workshop Project International Registry of HLA Epitopes

<http://www.epregistry.com.br>

Duquesnoy RJ, Marrari M, **da M Sousa LCD**, de M. Barroso JRP, de SE Aita KM, da Silva AS and **do Monte SJH** : 16th IHIW: a website for the antibody-defined HLA epitope registry. Int J. Immunogenetics, 40: 54-59, 2013

Databases for HLA-ABC, HLADRB1/3/4/5, HLA-DQA-B, HLA-DPA-B and MICA

The HLA Epitope Registry has now
the updated class I eplet repertoire

Recording Antibody-Defined Epitopes in Registry

- Each database in the HLA Epitope Registry comprises a list of potential epitopes that can be predicted to serve as recognition sites for specific antibodies
- The intention is to document epitopes that have been verified experimentally with informative specific antibodies

Recording Antibody-Defined Epitopes in Registry

- Each database in the HLA Epitope Registry comprises a list of potential epitopes that can be predicted to serve as recognition sites for specific antibodies
- The intention is to document epitopes that have been verified experimentally with informative specific antibodies
- Depending on the completeness of antibody reactivity information, there are two categories of antibody-defined epitopes: “confirmed” or “provisional”
- Analogous to the original serological assignments of HLA antigens during the early international workshops, the provisional status will be upgraded to confirmed if sufficient additional verification becomes available

Recording Antibody-Defined Epitopes in Registry

- Each database in the HLA Epitope Registry comprises a list of potential epitopes that can be predicted to serve as recognition sites for specific antibodies
- The intention is to document epitopes that have been verified experimentally with informative specific antibodies
- Depending on the completeness of antibody reactivity information, there are two categories of antibody-defined epitopes: “confirmed” or “provisional”
- Analogous to the original serological assignments of HLA antigens during the early international workshops, the provisional status will be upgraded to confirmed if sufficient additional verification becomes available
- For each antibody-defined epitope a link on the webpage leads to detailed information about antibody reactivity patterns

Antibody Verification of HLA Epitopes

- Duquesnoy RJ, Marrari M, Mulder A, da M. Sousa LCD, Da Silva AS, do Monte SJH. 'First Report on the Antibody Verification of HLA-ABC Epitopes Recorded in the HLA Epitope Registry'. *Tissue Antigens* 83:391-400, 2014
 - 97 HLA-ABC antibody-verified epitopes have been recorded, 62 correspond to eplets and 35 are defined by eplets paired with other residue configurations
- Duquesnoy RJ, Marrari M, Tambur A, da Mata Sousa LCD, Da Silva AS, do Monte SJH. 'First Report on the Antibody Verification of HLA-DR, HLA-DQ and HLA-DP Epitopes Recorded in the HLA Epitope Registry' *Human Immunology* 75:1097-1103, 2014
 - 24 HLA-DRB1/3/4/5, 15 DQB, 3 DQA and 8 DPB antibody-verified epitopes
- Duquesnoy RJ, Marrari M, Mosteck J, Da Silva AS, da Mata Sousa LCD, do Monte SJH. 'First Report on the Antibody Verification of MICA Epitopes Recorded in the HLA Epitope Registry' *Int J Immunogenetics* 41: 370-377, 2014
 - 21 MICA epitopes verified with antibodies which have primarily been tested in Luminex assays with single alleles

Antibody-verified HLA-ABC Epitopes corresponding to Eplets and Recorded in the Registry (October 2014)

Antibody Eplet	Antibody Verified	Antibody Eplet	Antibody Verified	Antibody Eplet	Antibody Verified	Antibody Eplet	Antibody Verified	Antibody Eplet	Antibody Verified	Antibody Eplet	Antibody Verified
9D		62RN		69ATS		76ES		102HV		151AHE	
9F		62RNQ		69RA		76ESI Confirmed		103L		151AHV	
9H		62RNR Provisional		69RT		76ESN Confirmed		103M ₂		151ARV	
9S		62RR Provisional		69TNT Provisional		76ET		105S		151H	
9T		62RRN		70HT		76VDT		107W Confirmed		152A	
9Y		62RTN		70IAQ Provisional		76VRN Provisional		109F		152E	
11AM		63EI		70QT		76VS		109FE		152HA	
11AV		63EK		71ATD Provisional		77NGT		113HD		152RA	
11SM		63EN		71ATN		77SRN		113HN		152RE	
16S		63ER		71ATS		79GT		113YD		151RR Provisional	
17RS Provisional		63ERN		71HS		80I Confirmed		113YH		152RT	
17WR		63NN		71KA		80K Provisional		113YN		152RW Provisional	
19K		65GK Confirmed		71QS		80N Confirmed		113YQ		152V	
1C		65QIA Confirmed		71SA Confirmed		80T		113YR		156DA Provisional	
21H Confirmed		65QKR		71STN		80TA		116D		156L	
24S		65QNR		71TTD		80TL		116F		156LA	
24T		65RA		71TTN		80TLR Provisional		116L		156QA	
30G		65RK		71TTS Provisional		81ALR		116S		156R	
32L		65RNA		72QTD		82LR Confirmed		116Y		156RA	
35Q		66I		73AN		90AR		11SV		156WA	
41T Confirmed		66IC		73AS		90D Provisional		127K Confirmed		158T Confirmed	
43RRM		66ICT		73ID		94II		131S		161D Confirmed	
44KM ₃ Confirmed		66IF		73TD		94IL		138K Provisional		162DLS	
44RG		66IS		73TDA		94TF		138MI		162GLS	
44RM		66IY		73TDE		94TI		142ITQ		163E	
44RMA Confirmed		66IYT		73TDV		94TL		143S		163EW Confirmed	
44RME		66K		73TN		94TV		144K Provisional		163L	
44RT Provisional		66KA		73TS		94TW		144KA		163LE	
45EE		66KAH		73TV		97I		144KR Confirmed		163LG	
45KE		66N		73TVD		97M		144QL Confirmed		163LW Confirmed	
56R Confirmed		66NAH		73TVN		97N		144TKH		163R Provisional	
62EE Provisional		66NAQ		73TVS Provisional		97R		145HT		163RG Confirmed	
62GE Confirmed		66NIM		73TY		97S		145KHA Provisional		163RW Confirmed	
62GK ₂ Confirmed		66NV		76ANT Confirmed		97T		145RT Provisional		163T	
62GRN Confirmed		66RKQ		76ED		97V		145L Confirmed		163TEW	
62LQ Confirmed		69AA Confirmed		76EG		97W		149AH		163TG	
62QE		69AQT		76EN		99F		149TAH Confirmed		166DG Confirmed	
62REN		69AT		76ENI		99S		149AAH Confirmed		166ES	
62RER		69ATD		76ENR		99Y		150AH		170RH	
62RK		69ATN		76ENT		102DV		151AHA Confirmed		170RY	

Antibody-verified HLA-ABC Epitopes corresponding to Eplet Pairs and Recorded in the Registry (October 2014)

HLA-A,B,C Epitope	Polymorphic Residue Description	HLA-A,B,C Epitope	Polymorphic Residue Description
44RT+69TNT	44R45T46E paired with 69T70N71T	82LR+138T	79R82L83R paired with 138T
62QE+56G	62Q63E65R66N paired with self56G	82LR+144Q	79R82L83R paired with 144Q
62QE+151H	62Q63E65R66N paired with self151H	82LR+144QR	79R82L83R paired with 144Q145R
62RN+163TW	62R63N paired with 163T166W	82LR+145R	79R82L83Rpaired with 145R
65QIA+76ESN	65Q66I69A paired with 76E77S80N	82LR+145RA	79R82L83R paired with145R149A
65QK+76VS	62R65Q66K67Y69R paired with self 76VS	90D+138M	90D91G paired with 138M
65RNA+ 80I	65R66N69A paired with self80I82L83R	131S+163LW	131S paired with 163L166W
66IF+163TW	66I67F69T paired with 163T166W	131S+163T	131S paired with 163T166W
69AA+ 80N	69A71A paired with self80N	138MI+79GT	138M142I143T145R paired with 79G80T
69TNT+80N	69T70N71T paired with 80N	143S+76ESN	142I143S144Q145R paired with self 76ESN
76ENI+62RE	76E77N79R80I paired with 62R62E	144K+76DVT	144K paired with 76DVT
79GT+19E	79G80T81L82R83G paired with 19E	144KR+151H	142I143T144K145R paired with self151H
80I+65QI	79R80I81A82L83R paired with 65Q66I	163EW+73TE	162G163E166E167W paired with 73T76E
80I+90A	79R80I81A82L83R paired with 90A or 145T	163LW+65QI	162G163L166E167Wpaired with self 65Q66I
80I+151RE	79R80I81A82L83R paired with 151R152E	163TEW+65QI	163T166E167W paired with 65Q66i
82LR+90A	79R82L83R paired with 90A	163TEW+103L	163T166E167Wpaired with 103L
82LR+138M	79R82L83R paired with 138M		

All have provisional status

Example of a record of an antibody-verified HLA-A epitope

	PANEL	OL MFI	GP MFI
	Positive control	9935	15973
	Negative control	31	127
Immunizing allele:	144TKH+ A*02:01	16165	17273
	144TKH+ A*02:02	nt	17579
	144TKH+ A*02:03	21200	17606
	144TKH+ A*02:05	nt	18235
	144TKH+ A*02:06	22505	nt
	144TKH+ A*68:01	17886	18358
	144TKH+ A*68:02	18034	17741
	144TKH+ A*69:01	19158	12210
	Self Alleles	35 ± 16	70 ± 21
	Other 144TKH-negative Alleles	31 ± 10	89 ± 30

Example of a record of an antibody-verified HLA-B epitope

Epitope: 41T (confirmed) on B13,B40,B41,B44,B45,B47,B49,B50

Report Date: 21Feb2012

Duquesnoy et al *Human Immunology* 73: 267-77, 2012

Pregnancy induced human monoclonal antibody ROU9A6 (IgG) from Arend Mulder (Leiden):


Antibody producer: HLA-A1, A25; B8, B18; Cw7.

Immunizer:HLA-A2; B44, B60; Cw5, w10.

El-Awar et al, *Human Immunology* 68: 170-180, 2007:

Serum eluate W7248.BO (absorbed with B*44:02 cells, no information on HLA types of antibody producer/immunizer)
Equivalent to TerEp #21

PANEL	OL MFI	GP MFI
Positive control	9488	nt
Negative control	9	nt
41T+ B*13:01	15174	nt
41T+ B*13:02	19404	18245
41T+ B*40:01	10610	22974
41T+ B*40:02	18522	21031
41T+ B*40:06	16626	nt
41T+ B*41:01	20733	21957
41T+ B*44:02	18295	7464
Immunizing allele: 41T+ B*44:03	18305	14695
41T+ B*45:01	21450	19061
41T+ B*47:01	17208	20014
41T+ B*49:01	13830	15192
41T+ B*50:01	17934	12216
Self Alleles	11 ± 2	68 ± 12
Other 41T-negative Alleles	76 ± 232	81 ± 41



Fernandez-Vina et al, in *Immunobiology of the Human MHC* vol. I, J.A. Hansen (Ed), pp 890-931, 2006:

ROU9A6 (13-WS 0120) positive with HLA-B*13:01, *13:02, *13:04, *40:01, *40:02, *40:04, *40:06, *40:11, *40:19, *40:30, *41:01, *41:02, *44:02, *44:03, *44:04, *44:21, *45:01, *47:01, *47:02, *47:03, *49:01, *50:01, *50:02 in CDC assays in 13th IHWS.

Example of a record of an antibody-verified HLA-C epitope

	PANEL	OL MFI	GP MFI
	Positive control	nt	nt
	Negative control	10	1
	219W+ C*01:02	13483	7444
	219W+ C*03:02	11733	nt
Immunizing allele:	219W+ C*03:03	13204	9740
	219W+ C*03:04	14135	10896
	219W+ C*04:01	7073	4549
	219W+ C*04:03	nt	7939
	219W+ C*14:02	12176	10036
	219W+ C*18:01	nt	10985
	219W+ C*18:02	12659	nt
	Self Alleles	83 ± 110	11 ± 1
	Other 219W-negative Alleles	109 ± 259	9 ± 3

Newly Antibody-Verified Epitopes

- For a successful registry we need contributions by HLA professionals who have identified antibodies specific for new and not so well described epitopes
- From the website one can download instructions how to submit information about epitope-specific antibody reactivity